

Max-550E Sensors

SENSING

ANALYSIS

DELIVERY

*"Same great sensor, but perfected
for specific applications"*



MaxBlend Family
Flow Diverter



MAX-550E FOR MAXO2 ME: P/N R140P02

- Works with all Max-250 Accessories
- Flow Diverter will not fit in MaxBlend Family

MAX-550E FOR MAXBLEND FAMILY P/N R140P02-001

- Works with the MaxBlend Family of Products
(MaxBlend 2, MaxBlend Lite, BlenderBuddy 2)
- Flow Diverter will not fit in standard Max-250 accessories

SO, WHY THE DIFFERENCE?

The **MaxBlend Family** of products are equipped with a bleed-control, which allows the user to provide **better accuracy** and **waste less gas**. Blenders require a minimum flow to be accurate. By adjusting the toggle to the appropriate position, the user will allow the gas to flow properly through the blender which **minimizes wasted gas to 0.1 LPM** when closed. Because of this, we had to fine-tune our Max-550E sensor to accurately sample the gas with such small traces of gas.

In Summary - *We needed to create a flow-diverter that can take an accurate reading with such a small bleed (less wasted gas).*

PART NUMBERS

FOR MAXBLEND FAMILY:

Max-550E, MaxBlend . . . R140P02-001
Flow Diverter, MaxBlend . R228P41-001

FOR ANALYZER FAMILY:

Max-550E, Analyzer R140P02
Flow Diverter, Analyzer . . . R110P10-001

ML# 0002 REV. A