

- **Compatible with the VDR-4 Respirator**

- **Features**

- FDA Registered; CE certified
- Lower audible noise level
- The volume regulator is designed to operate without a balloon
- Harness connections are specially designed to reduce risk and facilitate use
- Instructions for Use & Quick Reference Guide are included



- **Product**

Disposable, single patient use  
Non-sterile, ready-to-use from package  
Adult, pediatric, and neonate  
Duration of use: up to seven (7) days

- **Specifications**

Operates at frequencies of 60 – 900 cycles/min  
Accommodates a Flowrate up to 170L/min  
ISO 5356 Compliant connection port(s) 15mm ID & 22mm OD


- **Materials**

Medical grade plastic; compliant with biocompatibility requirements  
Latex free

- **Standard Packaging**

Units / Bag / Box 1 Unit / Bag / Box  
Unit / Case 5 Units / Case

*For further product information, please visit  
CircuitLife<sup>®</sup> website at [www.circuit-life.com](http://www.circuit-life.com)*



Provided By:  
  
TECHNOLOGY OF THE FUTURE, TODAY™  
866-624-3952  
[www.thinkMDI.net](http://www.thinkMDI.net)