



Making Safe Simple®

SĀFEN AWARDED THE 2023
"Weezie's Angel Healthcare Hero Award"
and "J. Pat O'Connell Technology
and Innovation Award"



LOUISE H. BATZ
PATIENT SAFETY FOUNDATION

GōTo® Tags



"I believe Sāfen® tags will reduce the incidence of medical errors."

– PK Carlton, Jr, Lt Gen (Ret) USAF, MD, FACS

Sāfen® GōTo® tags identify the best line for emergency medications.¹

- **Improve communication** and **situational awareness** for staff, patient and patient's family and support system
- Enter the room, or look from window and **immediately** know which IV line is for emergency use.

Intentionally designed to be a **simple, intuitive, "No Tech" solution** needing no alarms or computer interface.

Up to 53% of patients receiving IV medication will experience a medication error during a hospital stay.² Studies identify that up to 58% of these errors are related to incomplete or incorrect IV line labeling.³

1. In simulation-based testing, participants administered emergency medications quicker and more accurately when using the Sāfen GōTo™ Tag versus tape and marker.

2. Keers RN, Williams SD, Cooke J, Ashcroft DM. Prevalence and Nature of Medication Administration Errors in Health Care Settings: A Systematic Review of Direct Observational Evidence. *Annals of Pharmacotherapy* (2013);47(2):237-256.

3. Summa-Sorgini, Claudia, et al. "Errors associated with IV infusions in critical care." *The Canadian Journal of Hospital Pharmacy*. 65.1 (2012): 19.

SÄFEN AWARDED THE 2023
 "Weezy's Angel Healthcare Hero Award"
 and "J. Pat O'Connell Technology
 and Innovation Award"



LOUISE H. BATZ
 PATIENT SAFETY FOUNDATION

The Sāfen® Solution



GōTo® Tag Features:

- **Distinctive shape and color** for a quick visual of where the patient's emergency line is
- **Pre-printed text** on the front side eliminates handwriting and abbreviation confusion
- Arrow icons provide a **visual reminder** of the line and direction of fluid flow
- GōTo and arrows are raised, allowing the use of **touch/texture (haptics)** through gloved hands
- **Space to write** IV fluid specifics on the backside with a pen or marker

For Use in Multiple Patient Care Locations:



- Hospital (Emergency Department, Intensive Care Units, Medical-Surgical, Specialty Care Units, Operating Room)
- Prehospital
- Interfacility Transports
- Outpatient
- Home Health

GōTo® Tag Ordering Information:

Part # 88100 100 units/box

Learn more and/or request a demo:
info@safenmedical.com

FDA Registered Class I GMP Exempt tags are single use and disposable.
 Patented and Patents Pending. Made in the U.S.