

NEB BLOCK



The Neb Block makes your job easier.

It Increases Productivity ...

Neb treatments, resuscitation bags or equipment power are all available with this unique product. All equipment can be kept connected and ready for use when needed.

It Saves Time ...

You simply need to push the switch from “off” to “on” to deliver a neb treatment or to utilize a resuscitation bag.

It's Preset ...

The outlet barb is available in 6 pre-set flows for easy delivery of treatment.

It's Versatile ...









Need power for equipment? The Neb Block is equipped with a 50 psi power take off to allow use of any equipment that requires 50 psi power.

It's Convenient ...

The Neb Block is available with either chrome, compact or aluminum flowmeters. Connections for all wall outlet styles are available.

ORDERING INFORMATION (NEB BLOCKS)



Inlet Connection	Neb Drive	Neb Block	w/Chrome Body Flowmeter	w/Aluminum Flowmeter	w/Compact Flowmeter
No Adaptor (1/8 NPTF)	6-7 LPM	NB6060	1MFA6060	5MFA6060	8MFA6060
	7-8 LPM	NB6070	1MFA6070	5MFA6070	8MFA6070
	8-9 LPM	NB6080	1MFA6080	5MFA6080	8MFA6080
	10-11 LPM	NB6100	1MFA6100	5MFA6100	8MFA6100
	14-16 LPM	NB6150	1MFA6150	5MFA6150	8MFA6150
	24-26 LPM	NB6240	1MFA6240	5MFA6240	8MFA6240
 DISS Female Hex Nut	6-7 LPM	NB6062	1MFA6062	5MFA6062	8MFA6062
	7-8 LPM	NB6072	1MFA6072	5MFA6072	8MFA6072
	8-9 LPM	NB6082	1MFA6082	5MFA6082	8MFA6082
	10-11 LPM	NB6102	1MFA6102	5MFA6102	8MFA6102
	14-16 LPM	NB6152	1MFA6152	5MFA6152	8MFA6152
	24-26 LPM	NB6242	1MFA6242	5MFA6242	8MFA6242
 DISS Female Hand Tight	6-7 LPM	NB6063	1MFA6063	5MFA6063	8MFA6063
	7-8 LPM	NB6073	1MFA6073	5MFA6073	8MFA6073
	8-9 LPM	NB6083	1MFA6083	5MFA6083	8MFA6083
	10-11 LPM	NB6103	1MFA6103	5MFA6103	8MFA6103
	14-16 LPM	NB6153	1MFA6153	5MFA6153	8MFA6153
	24-26 LPM	NB6243	1MFA6243	5MFA6243	8MFA6243
 DISS Male	6-7 LPM	NB6064	1MFA6064	5MFA6064	8MFA6064
	7-8 LPM	NB6074	1MFA6074	5MFA6074	8MFA6074
	8-9 LPM	NB6084	1MFA6084	5MFA6084	8MFA6084
	10-11 LPM	NB6104	1MFA6104	5MFA6104	8MFA6104
	14-16 LPM	NB6154	1MFA6154	5MFA6154	8MFA6154
	24-26 LPM	NB6244	1MFA6244	5MFA6244	8MFA6244
 Ohmeda	6-7 LPM	NB6065	1MFA6065	5MFA6065	8MFA6065
	7-8 LPM	NB6075	1MFA6075	5MFA6075	8MFA6075
	8-9 LPM	NB6085	1MFA6085	5MFA6085	8MFA6085
	10-11 LPM	NB6105	1MFA6105	5MFA6105	8MFA6105
	14-16 LPM	NB6155	1MFA6155	5MFA6155	8MFA6155
	24-26 LPM	NB6245	1MFA6245	5MFA6245	8MFA6245
 Chemetron (NCG)	6-7 LPM	NB6066	1MFA6066	5MFA6066	8MFA6066
	7-8 LPM	NB6076	1MFA6076	5MFA6076	8MFA6076
	8-9 LPM	NB6086	1MFA6086	5MFA6086	8MFA6086
	10-11 LPM	NB6106	1MFA6106	5MFA6106	8MFA6106
	14-16 LPM	NB6156	1MFA6156	5MFA6156	8MFA6156
	24-26 LPM	NB6246	1MFA6246	5MFA6246	8MFA6246
 Oxequip	6-7 LPM	NB6067	1MFA6067	5MFA6067	8MFA6067
	7-8 LPM	NB6077	1MFA6077	5MFA6077	8MFA6077
	8-9 LPM	NB6087	1MFA6087	5MFA6087	8MFA6087
	10-11 LPM	NB6107	1MFA6107	5MFA6107	8MFA6107
	14-16 LPM	NB6157	1MFA6157	5MFA6157	8MFA6157
	24-26 LPM	NB6247	1MFA6247	5MFA6247	8MFA6247
 Puritan-Bennett	6-7 LPM	NB6068	1MFA6068	5MFA6068	8MFA6068
	7-8 LPM	NB6078	1MFA6078	5MFA6078	8MFA6078
	8-9 LPM	NB6088	1MFA6088	5MFA6088	8MFA6088
	10-11 LPM	NB6108	1MFA6108	5MFA6108	8MFA6108
	14-16 LPM	NB6158	1MFA6158	5MFA6158	8MFA6158
	24-26 LPM	NB6248	1MFA6248	5MFA6248	8MFA6248
 Schrader	6-7 LPM	NB6069	1MFA6069	5MFA6069	8MFA6069
	7-8 LPM	NB6079	1MFA6079	5MFA6079	8MFA6079
	8-9 LPM	NB6089	1MFA6089	5MFA6089	8MFA6089
	10-11 LPM	NB6109	1MFA6109	5MFA6109	8MFA6109
	14-16 LPM	NB6159	1MFA6159	5MFA6159	8MFA6159
	24-26 LPM	NB6249	1MFA6249	5MFA6249	8MFA6249

Custom assemblies available upon request.



Precision Medical, Inc.
 300 Held Drive • Northampton, PA 18067
 Phone: 610-262-6090 Toll-Free Phone: 1-800-272-7285
 Fax: 610-262-6080 Toll-Free Fax: 1-800-352-1240
 Visit us on the Web: www.precisionmedical.com