


# air-Q®

## Intubating Laryngeal Airways (ILA)

The everyday airway that's ready for the unexpected.

Striving  
to be  
**#1 IN PATIENT  
SATISFACTION**



Competitive ILA



air-Q with  
Epiglottis Elevator



Salter's expanding range of innovative products for airway management and capnography set the standard of performance and patient safety. Our products innovative features are designed by clinicians for clinicians for a difference you can see.




air-Q® Intubating Laryngeal Airways

**SALTER LABS®**  
DESIGN INNOVATIONS FOR  
RESPIRATORY AND AIRWAY CARE


# Perfect Placement Made Easy

Our unique intubating laryngeal airway includes an epiglottis elevator and raised heel to optimize airway placement and seal.



## air-Q Intubating SGA Airway provides:

- Anatomically correct placement of the masked airway
- Lifting of the epiglottis to maximize access
- Guide ramp to assist with ETT insertion and placement



# A Family of ILA's for All Your Needs

air-Q comes in 3 variants of intubating laryngeal airways that deliver safety and confidence in airway management, no matter what clinical challenges you face.

.....



## air-Q Standard

Traditional ILA for everyday case use that provides exceptional ease of ETT placement when unexpected intubation is required.



## air-Qsp (Self-Pressurizing)

Utilizes your positive pressure ventilatory source to automatically inflate the cuff to form a seal without manual inflation. Ideal for spontaneous trigger ventilation.



## air-Q Blocker

Traditional ILA features with added gastric access for cases where gastric reflux may be an issue. Ideal for use in the delivery room, emergency department and bariatric patients where ILA use is necessary.

### All air-Q products come standard with:

- Color coded cap to for quick sizing and is removable to accept a wide range of ETTs
- Exclusive Epiglottis Ramp helps to maximize visualization and minimize airway resistance
- ET tube ramp helps direct scopes and ETT's directly to the trachea

**air-Q Standard and air-Qsp available in single patient or reusable styles**

# Sizing Configurations, Kits and Accessories

air-Q is available in a wide range of sizes and kits to meet your needs for hospital and EMS use.



air-Q Disposable Intubating Laryngeal Airways (10/Box)

|                                                                                     | Size | Color  | Item    | Size/Description     | Maximum O.E.T.T. Size | Patient Weight |
|-------------------------------------------------------------------------------------|------|--------|---------|----------------------|-----------------------|----------------|
|    | 1.0  | Blue   | 10-3010 | air-Q 1.0 Disposable | 4.5                   | <7 kgs         |
|    | 1.5  | Green  | 10-3015 | air-Q 1.5 Disposable | 5.0                   | 7-17 kgs       |
|    | 2.0  | Orange | 10-3020 | air-Q 2.0 Disposable | 5.5                   | 17-30 kgs      |
|   | 2.5  | Yellow | 10-3025 | air-Q 2.5 Disposable | 6.5                   | 30-50 kgs      |
|  | 3.5  | Red    | 10-3035 | air-Q 3.5 Disposable | 7.5                   | 50-70 kgs      |
|  | 4.5  | Purple | 10-3045 | air-Q 4.5 Disposable | 8.5                   | 70-100 kgs     |

air-Q Reusable Intubating Laryngeal Airways (1/Pkg)

|                                                                                     | Size | Color  | Item    | Size/Description   | Maximum O.E.T.T. Size | Patient Weight |
|-------------------------------------------------------------------------------------|------|--------|---------|--------------------|-----------------------|----------------|
|  | 0.5  | Pink   | 10-1001 | air-Q 0.5 Reusable | 3.5                   | <4 kgs         |
|  | 1.0  | Blue   | 10-1010 | air-Q 1.0 Reusable | 4.5                   | 4-7 kgs        |
|  | 1.5  | Green  | 10-1015 | air-Q 1.5 Reusable | 5.0                   | 7-17 kgs       |
|  | 2.0  | Orange | 10-1020 | air-Q 2.0 Reusable | 5.5                   | 17-30 kgs      |
|  | 2.5  | Yellow | 10-1025 | air-Q 2.5 Reusable | 6.5                   | 30-50 kgs      |
|  | 3.5  | Red    | 10-1035 | air-Q 3.5 Reusable | 7.5                   | 50-70 kgs      |
|  | 4.5  | Purple | 10-1045 | air-Q 4.5 Reusable | 8.5                   | 70-100 kgs     |

# Sizing Configurations, Kits and Accessories



air-Q Disposable Airways/Lubricant/Syringe (10/Box)

|                                                                                   | Size | Color  | Item    | Size/Description                                | Maximum O.E.T.T. Size | Patient Weight |
|-----------------------------------------------------------------------------------|------|--------|---------|-------------------------------------------------|-----------------------|----------------|
|  | 1.5  | Green  | 10-3016 | 1.5 air-Q Disposable w/Lubricant & 10cc Syringe | 4.5                   | 7-17 kgs       |
|  | 2.0  | Orange | 10-3021 | 2.0 air-Q Disposable w/Lubricant & 10cc Syringe | 5.5                   | 17-30 kgs      |
|  | 2.5  | Yellow | 10-3026 | 2.5 air-Q Disposable w/Lubricant & 10cc Syringe | 6.5                   | 30-50 kgs      |
|  | 3.5  | Red    | 10-3036 | 3.5 air-Q Disposable w/Lubricant & 10cc Syringe | 7.5                   | 50-70 kgs      |
|  | 4.5  | Purple | 10-3046 | 4.5 air-Q Disposable w/Lubricant & 10cc Syringe | 8.5                   | 70-100 kgs     |

air-Qsp Disposable Intubating Laryngeal Airways (10/Box)

|                                                                                     | Size | Color  | Item    | Size/Description                                                       | Maximum O.E.T.T. Size | Patient Weight |
|-------------------------------------------------------------------------------------|------|--------|---------|------------------------------------------------------------------------|-----------------------|----------------|
|  | 1.0  | Blue   | 10-4010 | 1.0 air-Qsp (Self-Pressurizing) Disposable Intubating Laryngeal Airway | 4.5                   | 4-7 kgs        |
|  | 1.5  | Green  | 10-4015 | 1.5 air-Qsp (Self-Pressurizing) Disposable Intubating Laryngeal Airway | 5.0                   | 7-17 kgs       |
|  | 2.0  | Orange | 10-4020 | 2.0 air-Qsp (Self-Pressurizing) Disposable Intubating Laryngeal Airway | 5.5                   | 17-30 kgs      |
|  | 2.5  | Yellow | 10-4025 | 2.5 air-Qsp (Self-Pressurizing) Disposable Intubating Laryngeal Airway | 6.5                   | 30-50 kgs      |
|  | 3.5  | Red    | 10-4035 | 3.5 air-Qsp (Self-Pressurizing) Disposable Intubating Laryngeal Airway | 7.5                   | 50-70 kgs      |
|  | 4.5  | Purple | 10-4045 | 4.5 air-Qsp (Self-Pressurizing) Disposable Intubating Laryngeal Airway | 8.5                   | 70-100 kgs     |


# Sizing Configurations, Kits and Accessories



| air-Qsp Reusable Intubating Laryngeal Airways (1/Pkg)                              |      |        |         |                                                                      |                       |                |
|------------------------------------------------------------------------------------|------|--------|---------|----------------------------------------------------------------------|-----------------------|----------------|
|                                                                                    | Size | Color  | Item    | Size/Description                                                     | Maximum O.E.T.T. Size | Patient Weight |
|   | 0.5  | Pink   | 10-6005 | 0.5 air-Qsp (Self-Pressurizing) Reusable Intubating Laryngeal Airway | 4.0                   | <4 kgs         |
|   | 1.0  | Blue   | 10-6010 | 1.0 air-Qsp (Self-Pressurizing) Reusable Intubating Laryngeal Airway | 4.5                   | 4-7 kgs        |
|   | 1.5  | Green  | 10-6015 | 1.5 air-Qsp (Self-Pressurizing) Reusable Intubating Laryngeal Airway | 5.0                   | 7-17 kgs       |
|   | 2.0  | Orange | 10-6020 | 2.0 air-Qsp (Self-Pressurizing) Reusable Intubating Laryngeal Airway | 5.5                   | 17-30 kgs      |
|   | 2.5  | Yellow | 10-6025 | 2.5 air-Qsp (Self-Pressurizing) Reusable Intubating Laryngeal Airway | 6.5                   | 30-50 kgs      |
|   | 3.5  | Red    | 10-6035 | 3.5 air-Qsp (Self-Pressurizing) Reusable Intubating Laryngeal Airway | 7.5                   | 50-70 kgs      |
|  | 4.5  | Purple | 10-6045 | 4.5 air-Qsp (Self-Pressurizing) Reusable Intubating Laryngeal Airway | 8.5                   | 70-100 kgs     |

| air-Q Blocker Intubating Laryngeal Airways Disposable (10/Box)                      |      |        |         |                                     |                       |
|-------------------------------------------------------------------------------------|------|--------|---------|-------------------------------------|-----------------------|
|                                                                                     | Size | Color  | Item    | Size/Description                    | Maximum O.E.T.T. Size |
|  | 2.5  | Yellow | 10-5025 | 2.5 air-Q Blocker Disposable Airway | 6.5                   |
|  | 3.5  | Red    | 10-5035 | 3.5 air-Q Blocker Disposable Airway | 7.5                   |
|  | 4.5  | Purple | 10-5045 | 4.5 air-Q Blocker Disposable Airway | 8.5                   |

# Sizing Configurations, Kits and Accessories




air-Q Blocker Disposable Airways/Lubricant/Syringe (10/box)

|                                                                                   | Size | Color  | Item    | Size/Description                                                                      | Maximum O.E.T.T. Size |
|-----------------------------------------------------------------------------------|------|--------|---------|---------------------------------------------------------------------------------------|-----------------------|
|  | 2.5  | Yellow | 10-5026 | 2.5 air-Q Blocker Disposable Airway Kit w/Lubricant & 10cc Syringe with Swivel Elbow. | 6.5                   |
|  | 3.5  | Red    | 10-5036 | 3.5 air-Q Blocker Disposable Airway Kit w/Lubricant & 10cc Syringe with Swivel Elbow. | 7.5                   |
|  | 4.5  | Purple | 10-5046 | 4.5 air-Q Blocker Disposable Airway Kit w/Lubricant & 10cc Syringe with Swivel Elbow. | 8.5                   |

air-Q Blocker Tube (10/box)

|                                                                                     | Size | Color  | Item    | Size/Description       | Use With                                      |
|-------------------------------------------------------------------------------------|------|--------|---------|------------------------|-----------------------------------------------|
|   | 2.5  | Yellow | 10-5027 | 2.5 air-Q Blocker Tube | 2.5 air-Q Blocker Disposable Airway - 10-5025 |
|  | 3.5  | Red    | 10-5037 | 3.5 air-Q Blocker Tube | 3.5 air-Q Blocker Disposable Airway - 10-5035 |
|  | 4.5  | Purple | 10-5047 | 4.5 air-Q Blocker Tube | 4.5 air-Q Blocker Disposable Airway - 10-5045 |



air-Q® Stylet Guide (10/box)

| air-Q Size | Connector Color | Item    | Max O.E.T.T. | Patient Weight | Stylet Size |
|------------|-----------------|---------|--------------|----------------|-------------|
| 1.0        | Blue            | 10-1006 | 4.5          | 4-7 kgs        | 00 (Small)  |
| 1.5        | Green           | 10-1006 | 5.0          | 7-17 kgs       | 00 (Small)  |
| 2.0        | Orange          | 10-1005 | 5.5          | 17-30 kgs      | 0 (Medium)  |
| 2.5        | Yellow          | 10-1005 | 6.5          | 30-50 kgs      | 0 (Medium)  |
| 3.5        | Red             | 10-1004 | 7.5          | 50-70 kgs      | 1 (Large)   |
| 4.5        | Purple          | 10-1004 | 8.5          | 70-100 kgs     | 1 (Large)   |



# Salter Airway Solutions

Salter's expanding range of innovative products for airway management and capnography set the standard of performance and patient safety. Our products' innovative features were designed by actual clinicians, from the first true advance in laryngoscope lighting in decades to flexible endotracheal tube tips and now introducing the air-Q® family of Intubating Laryngeal Airways.


Improve intubation safety with Salter's innovative airway management products



air-Q Masked  
Laryngeal Airway



IntuBrite  
Laryngoscopes




Flex-Tip  
Endotracheal Tubes


## SalterSTAT Capnography



Orange Series for Reflective  
Connector Systems



Blue Series for Mask  
Applications



Green Series for Luer-Style  
Connector Systems

1. Ebert TJ, Novallija J, Uhrich TD, Barney JA. The effectiveness of oxygen delivery and reliability of carbon dioxide waveforms: a crossover comparison of 4 nasal cannulae. *Anesth Analg*. 2015;120(2):342-438.
2. Data on file. Salter Labs.
3. Salter Labs Technical Report 006

Provided By:  
**Medical Dynamics**  
TECHNOLOGY OF THE FUTURE, TODAY™  
866-624-3952  
[www.thinkMDI.net](http://www.thinkMDI.net)

International Support  
Toll Free: 303-799-1990  
Email: [international@salterlabs.com](mailto:international@salterlabs.com)

**SALTER LABS®**  
DESIGN INNOVATIONS FOR  
RESPIRATORY AND AIRWAY CARE