


ENPLUG



ENPLUGs are designed for emergency use to keep stomas from closing when a gastrostomy device has fallen out. The ENPLUG packs are available in 2 different lengths and various French sizes in case the stoma has already started to close and the current French gauge tube no longer fits the stoma.

Features & Benefits

- Stops stoma sites from closing
- Saves the stoma until a new gastrostomy tube can be placed
- Includes a range of French gauge sizes in each pack to provide options even when the stoma has started to close
- Effective, safe and inexpensive
- Saves staff time in emergency situations
- Promptly dealing with gastrostomy emergencies prevents use of other incorrect or inappropriate feeding devices
- Made from high grade medical silicone for flexibility
- Tapered tip for ease of insertion



Direct Codes	Specification	UOS
EN-PLUG-S	ENPLUG Silicone Stoma Plug Set, 4cm, 10Fr, 12Fr, 14Fr & 16Fr	100 units
EN-PLUG-L	ENPLUG Silicone Stoma Plug Set, 7cm, 10Fr, 12Fr, 14Fr & 16Fr	100 units
EN-PLUG-XL	ENPLUG Silicone Stoma Plug Set, 7cm, 18Fr & 20Fr	100 units