


Enteral pH



Features & Benefits

- Shelf life of 4 years
- Increments of 0.5
- Clear graduations between pH steps, especially between 4.5-5.0 & 5.0-5.5
- pH range of 2-9
- Can be stored from 4°C to 30°C



Available in packs of 25 for single patient use

Direct Codes	Specification	UOS
ENT-PHSPPX	Enteral pH Colour-Fixed Indicator Strips, 2.0-9.0	25 strips
ENT-PHX	Enteral pH Colour-Fixed Indicator Strips, 2.0-9.0	100 strips