

**PRODUCT  
CATALOGUE  
2019**



**Nasogastric  
Feeding Tubes**

**Gastrostomy  
& Jejunostomy  
Feeding**

**Syringes**

**Accessories**

## Nasogastric Feeding

- |   |   |   |
|---|---|---|
| <b>03</b> Nutricare                     | <b>07</b> NG Tube Safety Pack & CareTip Grip            | <b>11</b> Bridle Pro<br><b>AMT Bridle<sup>†</sup></b> |
| <b>04</b> Nutricare Infant Feeding Tube | <b>08</b> Enteral Drainage Bags                         | <b>12-13</b> Enteral pH & PHX-Act                     |
| <b>05</b> Carefeed                      | <b>09</b> Enteral Extension Sets                        | <b>14</b> pHX   |
| <b>06</b> Enteral Connector Brush       | <b>10</b> Adult Bridle<br><b>AMT Bridle<sup>†</sup></b> |   |

## Gastrostomy/Jejunostomy Feeding

- |   |   |  |
|---|---|--|
| <b>15</b> Initial Placement Kit & Dilator Kit           | <b>20-21</b> MiniONE <sup>†</sup> Balloon Button  | <b>25</b> MiniONE <sup>†</sup> Button Feeding Sets                     |
| <b>16</b> Balloon G-Tube <sup>†</sup> G7 Pack           | <b>22</b> G-Pak                                   | <b>26-27</b> Low Profile G-Jet <sup>†</sup> (Transgastric Jejunostomy) |
| <b>17</b> Capsule Monarch <sup>†</sup> Gastrostomy Tube | <b>23</b> MiniONE <sup>†</sup> Non-Balloon Button | <b>28</b> ENPLUG <sup>*</sup>  |
| <b>18-19</b> MINI Balloon Button                        | <b>24</b> MINI Balloon Button Feeding Sets        |  |

## Enteral Syringes


- |                                       |                                |
|---------------------------------------|--------------------------------|
| <b>29</b> Enteral Syringe Driver      | <b>32</b> Oral Syringes        |
| <b>30</b> Enteral Single Use Syringes | <b>33</b> Catheter Tip Syringe |
| <b>31</b> Enteral Home Use Syringes   | <b>33</b> Colostrum Gold Kit   |

## Accessories

- |   |                                     |  |   |
|---|-------------------------------------|--|---|
| <b>34</b> ENFit <sup>®</sup> Bottle Adaptor | <b>35</b> Filter Needle             | <b>36</b> Fixation Plaster                     | <b>37</b> AMT Introducer <sup>†</sup>                     |
| <b>34</b> Enteral Medicine Straw            | <b>35</b> Y-Port Adaptor            | <b>36</b> Oral Bottle Adaptor                  | <b>37</b> AMT Stiffener Guidewire <sup>†</sup>            |
| <b>34</b> Medicine Bottle Adaptor           | <b>36</b> Button Bear               | <b>36</b> Enteral Care Case & Demo Case        | <b>37</b> AMT Stiffener Guidewire Compatible <sup>†</sup> |
| <b>35</b> Enteral Syringe Cap               | <b>36</b> Enteral Halo <sup>*</sup> | <b>36</b> ENFit <sup>®</sup> to Funnel Adaptor | <b>37</b> AMT Stoma Measuring Devices <sup>†</sup>        |
| <b>35</b> ENFit <sup>®</sup> Female Cap     | <b>36</b> Side Board                | <b>36</b> Funnel to ENFit <sup>®</sup> Adaptor | <b>37</b> AMT Clamp <sup>†</sup>                          |

## Long Term NG Tubes & NJ Tubes

Nutricare ISOSAF is a fully radiopaque polyurethane tube with clear numbered cm markings. It is fitted with the new ENFit (ISO 80369-3) connector. Each pack includes the NG tube with guidewire\*, fixation plaster, record of insertion sticker and safety literature\*\*. NJ tube clearly identifiable with blue cap & tube markings.



### Nutricare Long Term NG Tube Order Codes

Size	Length	Direct Codes	NHS SC Code	Packaging	Unit of Sale
6Fr	50cm	NC06.050ISOSAF	FWM2403	1pc/pouch	5pcs./box
	80cm	NC06.080ISOSAF	FWM2404	1pc/pouch	5pcs./box
8Fr	50cm	NC08.050ISOSAF	FWM2405	1pc/pouch	5pcs./box
	80cm	NC08.080ISOSAF	FWM2406	1pc/pouch	5pcs./box
	92cm	NC08.092ISOSAF	FWM2407	1pc/pouch	5pcs./box
	120cm	NC08.120ISOSAF	FWM2408	1pc/pouch	5pcs./box
10Fr	92cm	NC10.092ISOSAF	FWM2409	1pc/pouch	5pcs./box
12Fr	92cm	NC12.092ISOSAF	FWM2410	1pc/pouch	5pcs./box
14Fr	92cm	NC14.092ISOSAF-ROW	FWM2752	1pc/pouch	5pcs./box
16Fr	92cm	NC16.092ISOSAF-ROW	FWM2753	1pc/pouch	5pcs./box

### Product Features

- USE UP TO **90 DAYS**
- FULLY RADIOPAQUE**
- NUMBERED EVERY CM**
- MALE ENFit CONNECTOR** (ISO-Saf)
- FIXATION PLASTER**
- LATEX FREE**
- NJ CLEARLY IDENTIFIED**


### Nutricare Long Term NJ Tube with Guidewire Order Codes

8Fr	140cm	NC08.140ISOSAF	FWM3177	1pc/pouch	5pcs./box
10Fr	140cm	NC10.140ISOSAF	FWM3178	1pc/pouch	5pcs./box

### Which Accessories & Related Products Can I Use?

- ENTERAL CONNECTOR BRUSH** (See page 06)
- ENTERAL DRAINAGE BAGS** (ISO-Saf) (See page 08)
- ENTERAL EXTENSION SETS** (ISO-Saf) (See page 09)
- Y-PORT ADAPTOR** (ISO-Saf) (See page 35)
- ENFit FEMALE CAP** (ISO-Saf) (See page 35)

**NUTRICARE  
INFANT  
CLEAR-**  
ISO 



### Infant Feeding Tube

Nutricare Infant Clear is a polyurethane tube with clear numbered cm marking. It is fitted with the new ENFit (ISO 80369-3) connector. Each single blister includes the NG tube and fixation plaster.

\*Blue stripe on tube is for illustrative purposes only

### Nutricare Infant Clear NG Tubes Order Codes

Size	Length	Direct Codes	NHS SC Code	Packaging	Unit of Sale
4Fr	50cm	NCI04.050ISOSAF	FTA331	1pc/pouch	25pcs./box
	75cm	NCI04.075ISOSAF	FTA332	1pc/pouch	25pcs./box
5Fr	50cm	NCI05.050ISOSAF	FTA333	1pc/pouch	25pcs./box
	75cm	NCI05.075ISOSAF	FTA334	1pc/pouch	25pcs./box
	125cm	NCI05.125ISOSAF	FTA335	1pc/pouch	25pcs./box
6Fr	50cm	NCI06.050ISOSAF	FTA336	1pc/pouch	25pcs./box
	75cm	NCI06.075ISOSAF	FTA337	1pc/pouch	25pcs./box
	125cm	NCI06.125ISOSAF	FTA338	1pc/pouch	25pcs./box
8Fr	50cm	NCI08.050ISOSAF	FTA339	1pc/pouch	25pcs./box
	75cm	NCI08.075ISOSAF	FTA340	1pc/pouch	25pcs./box
10Fr	75cm	NCI10.075ISOSAF	FTA341	1pc/pouch	25pcs./box

### Product Features

- RADIOPAQUE STRIPE**
- MALE ENFit CONNECTOR**  
ISO 
- FIXATION PLASTER**
- USE UP TO 90 DAYS**
- NUMBERED EVERY CM**
- LATEX FREE**

### Priming Volumes

Size	Length	Priming Volume
4Fr	50cm	0.25ml
	75cm	0.5ml
5Fr	50 cm	0.5ml
	75 cm	0.75ml
6Fr	50 cm	0.70ml
	75 cm	1.1ml
8Fr	50 cm	1.25ml
	75 cm	1.85ml



### Which Accessories & Related Products Can I Use?


**ENTERAL CONNECTOR BRUSH**  
  
See page 06

**ENTERAL DRAINAGE BAGS**  
  
See page 08

**ENTERAL EXTENSION SETS**  
  
See page 09


**Y-PORT ADAPTOR**  
  
See page 35

**ENFit FEMALE CAP**  
  
See page 35



## Short Term Feeding & Drainage Tube

Carefeed ISOSAF is made from PVC, has clear numbered cm markings and is fitted with the new ENFit (ISO 80369-3) connector. Each single blister includes the NG tube, fixation plaster & record of insertion sticker.



**CAREFEED**  
ISO Saf

### Carefeed Feeding & Drainage Tube Order Codes

Size	Length	Direct Codes	NHS SC Code	Packaging	Unit of Sale
4Fr	50cm	CF01.04.050ISOSAF	FWM2644	1pc/pouch	25pcs./box
	75cm	CF01.04.075ISOSAF	FWM2645	1pc/pouch	25pcs./box
5Fr	50cm	CF01.05.050ISOSAF	FWM2421	1pc/pouch	25pcs./box
	75cm	CF01.05.075ISOSAF	FWM2646	1pc/pouch	25pcs./box
6Fr	50cm	CF01.06.050ISOSAF	FWM2420	1pc/pouch	25pcs./box
	75cm	CF01.06.075ISOSAF	FWM2422	1pc/pouch	25pcs./box
8Fr	50cm	CF01.08.050ISOSAF	FWM2647	1pc/pouch	25pcs./box
	75cm	CF01.08.075ISOSAF	FWM2423	1pc/pouch	25pcs./box
	92cm	CF01.08.092ISOSAF	FWM2648	1pc/pouch	25pcs./box
10Fr	75cm	CF01.10.075ISOSAF	FWM2652	1pc/pouch	25pcs./box
	92cm	CF01.10.092ISOSAF	FWM2649	1pc/pouch	25pcs./box
12Fr	75cm	CF01.12.075ISOSAF	FWM2643	1pc/pouch	25pcs./box
	92cm	CF01.12.092ISOSAF	FWM2650	1pc/pouch	25pcs./box
14Fr	92cm	CF01.14.092ISOSAF	FWM2739	1pc/pouch	25pcs./box

### Product Features

- USE FOR 7 TO 10 DAYS\*
- FULLY RADIOPAQUE
- NUMBERED EVERY CM
- MALE ENFit CONNECTOR ISO Saf
- DEHP FREE
- LATEX FREE
- FIXATION PLASTER

\* Up to 7 days for Neonates

### Priming Volumes

Size	Length	Priming Volume	Size	Length	Priming Volume	Size	Length	Priming Volume
5Fr	50 cm	0.45ml	6Fr	50 cm	0.65ml	8Fr	50 cm	1.0ml
	75 cm	0.75ml		75 cm	1.0ml		75 cm	1.75ml

### Which Accessories & Related Products Can I Use?

ENTERAL DRAINAGE BAGS ISO Saf

See page 08

ENTERAL EXTENSION SETS ISO Saf

See page 09

Y-PORT ADAPTOR ISO Saf

See page 35

ENFit FEMALE CAP ISO Saf

See page 35



**ENTERAL  
CONNECTOR  
BRUSH**

**Enteral Connector Brush**

For keeping ENFit® connectors free of debris.



Connector Brush Order Code

Direct Code	NHS SC Code	Qty	Unit of Sale
ECB-NS	FWM3176	50	1

Which Accessories & Related Products Can I Use?

**NUTRICARE**  
ISO-Saf

See page 03

**NUTRICARE  
INFANT  
CLEAR-**  
ISO-Saf

See page 04

**CAREFEED**  
ISO-Saf

See page 05



## Award Winning NG Safety Insertion Material

INITIAL  
PLACEMENT  
NG TUBE  
SAFETY  
PACK

### NG Tube Safety Pack Order Codes

Direct Code	NHS SC Code
NGSP	FWM1827

### Product Features

NEX  
TAPE  
MEASURE

RECORD  
OF INSERTION  
STICKER

CareTip  
GRIP

\*Sold Separately

### Which Accessories & Related Products Can I Use?



## CareTip Grip Skin Fixation System

Unique skin fixation system designed specifically to hold tubing in place. e.g. NG, Feeding extensions or Jejunal tubes.

CareTip  
GRIP

### CareTip Grip Order Codes

	Direct Code	NHS SC Code	Unit of Sale
25 pairs	CTG-50-NS	FWM3200	50
100 in dispenser	CTG-100-NS	FWM3201	100

### Which Accessories & Related Products Can I Use?



**ENTERAL DRAINAGE BAGS**  
ISO 5Saf

**Enteral Drainage Bags**

Enteral Drainage Bags to fit any ENFit compliant tubes used for gastric drainage.

The 500ml and 1000ml bags include integral tubing. The bags are made from soft EVA material.



Male ENFit Connector



Funnel Connector

Enteral Drainage Bags Order Codes

With Male ENFit Connector To Distal End

Size	Codes	NHS SC Code	Packaging	Unit of Sale
50ml	DBISO-50	FWM2416	1pc/pouch	20pcs./box
100ml	DBISO-100	FWM2868	1pc/pouch	20pcs./box
250ml	DBISO-250	FWM2417	1pc/pouch	20pcs./box
500ml	DBISON-500	FWM2869	1pc/pouch	20pcs./box

With 50cm Inlet Tubing & Funnel Connector

Size	Codes	NHS SC Code	Packaging	Unit of Sale
500ml	DBISO-500	FWM2418	1pc/pouch	20pcs./box
1000ml	DBISO-1000	FWM2419	1pc/pouch	20pcs./box

Product Features

- CLAMP
- TETHERED CAP
- DEHP FREE
- ENFit COMPLIANT
- LATEX FREE

Which Accessories & Related Products Can I Use?

**NUTRICARE**  
ISO 5Saf

See page 03

**BALLOON G-TUBE**  
ISO 5Saf

See page 16

**G-JET**  
LOW PROFILE  
EXTENDING  
ENTERAL TUBE

See page 26-27

**ENTERAL**  
ISO 5Saf

See page 30


**ENTERAL**  
ISO 5Saf | Home

See page 31



## Enteral Extension Sets

Extension set/syringe driver connection set, male ENFit connector (ISO 80369-3) and female ENFit connector.



ENTERAL  
EXTENSION  
SETS  
ISO 80369-3

### Enteral Extension Sets Order Codes

Size	Direct Codes	NHS SC Code	Packaging	Unit of Sale
75cm	ESCF.075-ISOSAF	FWM2634	1pc/pouch	25pcs./box
150cm	ESCF.150-ISOSAF	FWM2635	1pc/pouch	25pcs./box

### Product Features



### Priming Volumes

Length	Priming Volume
75cm	1.3ml
150cm	2.6ml

### Which Accessories & Related Products Can I Use?





## Adult Bridle

When NG tubes won't stay in place, Bridle affords clinicians valuable time to determine the most appropriate nutritional program for their patients. Easily placed by passing tiny catheter-mounted magnets, Bridle's clip grips the NG Tube to positively secure it in place. Early and potentially unnecessary PEG tube placements can be avoided. Available in 8-18Fr sizes.

**AMT BRIDLE™**  
72% reduction  
in NG Tube  
pull-outs...



*Preserving "precious" naso-jejunal tubes...  
and improving patient outcomes.*

### Standard Bridle Order Codes

Size	Clip Colour	Direct Code	NHS SC Code	Unit of Sale
8Fr	Clip is WHITE	4-4108	FWM993	1
10Fr	Clip is TEAL	4-4110	FWM1120	1
12Fr	Clip is BLUE	4-4112	FWM995	1
14Fr	Clip is Lt. YELLOW	4-4114	FWM1533	1
16Fr	Clip is Lt. YELLOW	4-4116	FWM1388	1
18Fr	Clip is BLUE	4-4118	FWM2182	1

### Kit Contents

- Retrieval probe
- Bridle catheter with attached umbilical tape & stylet guide
- Attached feeding tube retaining clip

### Spare Clip Order Codes

	Direct Code	NHS SC Code	Unit of Sale
8Fr	4-4508	FWM1565	5
10Fr	4-4510	FWM1566	5
12Fr	4-4512	FWM1567	5
14Fr	4-4514	FWM2183	5
16Fr	4-4516	FWM1568	5

## Bridle Pro

Introducing the new AMT Bridle Pro™. AMT Bridle Pro is supplied with pre attached Pro Clip or pre attached soft over moulded pro range clips. AMT Bridle Pro tubing uses advanced polymer materials creating a smaller surface area for ease of placement and easier cleaning. Suitable for adults and paediatric.

- Available for all tubes 5Fr - 18Fr
- AMT Bridle Pro has a 60% less surface area resulting in greater comfort.



“

*The AMT Bridle™ is a useful tool in paediatric enteral nutrition and we would recommend it to other centres*”


”

The Great North Children's Hospital



### Bridle Pro Order Codes

Size	Product Description	Direct Code	NHS Code	Unit of Sale
5Fr-6Fr	PURPLE - Pro Range Kit <sup>1</sup>	4-420506	FWM3276	1 per box
8Fr-10Fr	TEAL - Pro Range Kit <sup>1</sup>	4-420810	FWM3277	1 per box
8Fr	WHITE - Pro Kit	4-4208	FWM3164	1 per box
10Fr	TEAL - Pro Kit	4-4210	FWM3165	1 per box
12Fr	BLUE - Pro Kit	4-4212	FWM3166	1 per box
14Fr	YELLOW - Pro Kit	4-4214	FWM3167	1 per box
16Fr	YELLOW - Pro Kit	4-4216	FWM3168	1 per box
18Fr	LIGHT BLUE - Pro Kit	4-4218	FWM3169	1 per box



Pro Range Clip  
8Fr-10Fr



Pro Clip  
10Fr

### Spare Clip Order Codes

5Fr-6Fr	PURPLE - Pro Spare Clip	4-460506	FWM3279	5 per box
8Fr-10Fr	TEAL - Pro Spare Clip	4-460810	FWM3275	5 per box
8Fr	WHITE - Pro Spare Clip	4-4608	FWM3268	5 per box
10Fr	TEAL - Pro Spare Clip	4-4610	FWM3269	5 per box
12Fr	BLUE - Pro Spare Clip	4-4612	FWM3270	5 per box
14Fr	YELLOW - Pro Spare Clip	4-4614	FWM3271	5 per box
16Fr	YELLOW - Pro Spare Clip	4-4616	FWM3272	5 per box
18Fr	LIGHT BLUE - Pro Spare Clip	4-4618	FWM3273	5 per box



5Fr-6Fr



8Fr-10Fr



12Fr-14Fr



16Fr-18Fr

<sup>1</sup> The Pro Range Clip's flexible over moulded inner section allows it to encompass multiple sizes into a single clip. Pro Range Clips 5Fr-6Fr and 8Fr-10Fr are smaller for paediatric patient use. \*Additional Sizes Coming Soon.





## How To Use Enteral pH & Troubleshooting

[www.gbukentral.com/videos](http://www.gbukentral.com/videos)



### Enteral-pH 2-9 Order Codes

Qty	Direct Code	NHS SC Code	Qty	Unit of Sales
25	ENT-PHSPPX	FWM3214	25	1
100	ENT-PHX	FWM3213	100	1

### Product Features

- CE FOR HUMAN GASTRIC ASPIRATE
- 0.5 INCREMENTS
- 4 YEAR SHELF LIFE
- SMALL SINGLE PATIENT PACK AVAILABLE

### Which Accessories & Related Products Can I Use?

**NUTRICARE**  
ISO 5834

See page 03

**NUTRICARE**  
INFANT CLEAR-ISO 5834

See page 04

**BALLOON G-TUBE**  
ISO 5834

See page 16

**ENTERAL**  
ISO 5834

See page 30

**ENTERAL**  
ISO 5834 Home

See page 31

## Innovative pH Strip Reader

A pH Strip reader that accurately determines the colour of Enteral pH strips and provides a printed record to attach to patient notes. A portable, safe and easy to use analyser that removes subjectivity from reading pH, leading to improved patient safety.



pHX-act

### pH Strip Reader Order Codes

Direct Codes	NHS SC Code	Unit of Sale
PHX-ACT	FWM3216	1

### pHX-act Printer Paper

Direct Codes	NHS SC Code	Unit of Sale
pHX-act-pp5	FWM3253	5 rolls/box

**A portable, safe and easy to use analyser that removes subjectivity from reading pH, leading to improved patient safety.**

### Product Features

- CE FOR HUMAN GASTRIC ASPIRATE
- EASY TO USE
- 0.5 INCREMENTS
- PORTABLE
- ACCURATE TO pH +/-0.2
- SAFE & ACCURATE
- NPSA ALERT  
NPSA/2011/PSA002
- HAS THE BENEFIT OF A PERMANENT RECORD OF THE MEASUREMENT
- TO BE USED WITH ENT-pHX & ENT-pHSPPX



## Introducing pH-X

A simple and effortless method for ascertaining the pH value of solutions with a colour scale you are familiar with.



### pHX Order Codes

Size	Codes	NHS SC Code	Packaging	Unit of Sale
0-6	pHX0-6	FWM3241	100pc/pot	1
2-9	pHX2-9	FWM3215	100pc/pot	1

### Product Features

- CE FOR HUMAN GASTRIC ASPIRATE
- 0-6 & 2-9
- 0.5 INCREMENTS
- HIGH QUALITY pH PAPER
- ACCURATE TO pH +/-0-2
- 100 PER POT
- NPSA ALERT  
NPSA/2011/PSA002

## INITIAL PLACEMENT KIT



### Comprehensive Kit

- Guidewire
- Introducer Needle, 18 Gauge
- Safety Scalpel, #11
- 5 Dilators (5, 7, 12, 16, 18, 22Fr)
- (1) Guidewire Compatible Introducer Stiffener
- 5cc Syringe
- Guidewire Compatible Balloon Stoma Measuring Device
- (4) 4" x 4" Gauze

### Comprehensive Kit Order Codes

Direct Code	NHS SC Code
IP-1218	FMW2177

Direct Code	NHS SC Code
IP-1218-S	FMW2178

Includes peel-away sheath on 18Fr dilator.

## DILATOR KIT



### Dilator Kit

- Guidewire
- (4) Dilators (7, 12, 16, 20Fr)
- Guidewire Compatible Introducer/Stiffener
- Introducer Needle, 18 Gauge


### Dilator Kit Order Codes

Direct Code	NHS SC Code
IP-DIL	FMW2179



**BALLOON  
G-TUBE**

ISO  Saf




**NEW AMT Balloon  
Gastrostomy Tube**

The AMT Balloon G-Tube is a Silicone replacement gastrostomy tube (G-Tube).

**G7 Pack**

1 x AMT Balloon Gastrostomy Tube and 1 x G-Pak.



Balloon G-Tube Order Codes

	Direct Code	Unit of Sale
12Fr	7-1210-I	1
14Fr	7-1410-I	1
16Fr	7-1610-I	1
18Fr	7-1810-I	1
20Fr	7-2010-I	1

Product Features



Which Accessories & Related Products Can I Use?




G7 Pack Order Codes

Direct Code	NHS SC Code	Contents
G7-7-1210-ISOSAF	FWM3121	1 x 7-1210-I + 1 X G-PAKISO
G7-7-1410-ISOSAF	FWM3122	1 x 7-1410-I + 1 X G-PAKISO
G7-7-1610-ISOSAF	FWM3123	1 x 7-1610-I + 1 X G-PAKISO
G7-7-1810-ISOSAF	FWM3124	1 x 7-1810-I + 1 X G-PAKISO








## Capsule Monarch Gastrostomy Tube

Unlike other G-Tubes, Capsule Monarch is not held in place by a balloon but uses an internal retention disc that is deployed when the Capsule Monarch® has been inserted into the stomach. Single ENFit ISO (80369-3) connector 12Fr & 14Fr. Y-Port Double ENFit connector 20Fr.



### Capsule Monarch Order Codes

	Direct Code	NHS SC Code	Unit of Sale
12Fr	9-C1210-I	FWM3125	1
14Fr	9-C1410-I	FWM3126	1
20Fr	9-C2010-I	FWM3127	1

#### Including Guidewire

12Fr	9-CG1210-I	FWM3175	1
14Fr	9-CG1410-I	FWM3174	1

### Product Features



### Capsule Monarch Accessories

	Capsule Monarch Accessories	Direct Code	NHS SC Code	Unit of Sale
12Fr	Monarch Replacement Single ENFit + Clamp	4-2012 ISOSAF	FWM2607	10
14Fr	Monarch Replacement Single ENFit + Clamp	4-2014 ISOSAF	FWM2608	10
20Fr	Monarch Replacement Single ENFit + Clamp	4-2020 ISOSAF	FWM2433	10



### Which Accessories & Related Products Can I Use?

**ENTERAL DRAINAGE BAGS**  
ISO Saf

See page 08

**ENTERAL EXTENSION SETS**  
ISO Saf

See page 09

**Y-PORT ADAPTOR**  
ISO Saf

See page 35

**ENFit FEMALE CAP**  
ISO Saf

See page 35

E5

**MINI**  
BALLOON BUTTON

### What is an E5?

A single MINI Balloon Button E5 with 5ml syringe.



### MINI Balloon Button E5 Order Codes

12Fr	Direct Code	NHS SC Code
0.8cm	E5-1208	FWM1550
1.0cm	E5-1210	FWM1551
1.2cm	E5-1212	FWM1552
1.5cm	E5-1215	FWM1553
1.7cm	E5-1217	FWM1554
2.0cm	E5-1220	FWM1555
2.3cm	E5-1223	FWM1556
2.5cm	E5-1225	FWM1557
2.7cm	E5-1227	FWM1558
3.0cm	E5-1230	FWM1559
3.5cm	E5-1235	FWM1560
4.0cm	E5-1240	FWM1592

14Fr	Direct Code	NHS SC Code
0.8cm	E5-1408	FWM1593
1.0cm	E5-1410	FWM1534
1.2cm	E5-1412	FWM1535
1.5cm	E5-1415	FWM1536
1.7cm	E5-1417	FWM1537
2.0cm	E5-1420	FWM1538
2.3cm	E5-1423	FWM1539
2.5cm	E5-1425	FWM1540
2.7cm	E5-1427	FWM1541
3.0cm	E5-1430	FWM1542
3.5cm	E5-1435	FWM1543
4.0cm	E5-1440	FWM1544
4.4cm	E5-1444	FWM1545
5.0cm	E5-1450	FWM1546
5.5cm	E5-1455	FWM1594
6.0cm	E5-1460	FWM1595
6.5cm	E5-1465	FWM1596

16Fr	Direct Code	NHS SC Code
1.0cm	E5-1610	FWM1597
1.2cm	E5-1612	FWM1598
1.5cm	E5-1615	FWM1599
1.7cm	E5-1617	FWM1600
2.0cm	E5-1620	FWM1601
2.3cm	E5-1623	FWM1602
2.5cm	E5-1625	FWM1603
2.7cm	E5-1627	FWM1604
3.0cm	E5-1630	FWM1605
3.5cm	E5-1635	FWM1606
4.0cm	E5-1640	FWM1608
4.4cm	E5-1644	FWM1607

18Fr	Direct Code	Availability
1.0cm	E5-1810	4-6 Weeks
1.2cm	E5-1812	4-6 Weeks
1.5cm	E5-1815	4-6 Weeks
1.7cm	E5-1817	4-6 Weeks
2.0cm	E5-1820	4-6 Weeks
2.3cm	E5-1823	4-6 Weeks
2.5cm	E5-1825	4-6 Weeks
2.7cm	E5-1827	4-6 Weeks
3.0cm	E5-1830	4-6 Weeks
3.5cm	E5-1835	4-6 Weeks
4.0cm	E5-1840	4-6 Weeks
4.4cm	E5-1844	4-6 Weeks

20Fr	Direct Code	Availability
1.0cm	E5-2010	4-6 Weeks
1.2cm	E5-2012	4-6 Weeks
1.5cm	E5-2015	4-6 Weeks
1.7cm	E5-2018	4-6 Weeks
2.0cm	E5-2020	4-6 Weeks
2.3cm	E5-2023	4-6 Weeks
2.5cm	E5-2025	4-6 Weeks
2.7cm	E5-2027	4-6 Weeks
3.0cm	E5-2030	4-6 Weeks
3.5cm	E5-2035	4-6 Weeks
4.0cm	E5-2040	4-6 Weeks

24Fr	Direct Code	Availability
1.5cm	E5-2415	4-6 Weeks
1.7cm	E5-2417	4-6 Weeks
2.0cm	E5-2420	4-6 Weeks
2.3cm	E5-2423	4-6 Weeks
2.5cm	E5-2425	4-6 Weeks
2.7cm	E5-2427	4-6 Weeks
3.0cm	E5-2430	4-6 Weeks
3.5cm	E5-2435	4-6 Weeks
4.0cm	E5-2440	4-6 Weeks

**MINI Balloon Button Feeding Sets**  
See page 24

**Non-stock item:** orders placed on request available in 4-6 weeks.

## MINI Balloon Button E7 Order Codes

12Fr	Direct Code	NHS SC Code
0.8cm	E7-1208 ISOSAF	FWM2567
1.0cm	E7-1210 ISOSAF	FWM2577
1.2cm	E7-1212 ISOSAF	FWM2588
1.5cm	E7-1215 ISOSAF	FWM2414
1.7cm	E7-1217 ISOSAF	FWM2413
2.0cm	E7-1220 ISOSAF	FWM2415
2.3cm	E7-1223 ISOSAF	FWM2626
2.5cm	E7-1225 ISOSAF	FWM2620
2.7cm	E7-1227 ISOSAF	FWM2621
3.0cm	E7-1230 ISOSAF	FWM2631
3.5cm	E7-1235 ISOSAF	FWM2632
4.0cm	E7-1240 ISOSAF	FWM2614

14Fr	Direct Code	NHS SC Code
0.8cm	E7-1408 ISOSAF	N/A
1.0cm	E7-1410 ISOSAF	FWM2578
1.2cm	E7-1412 ISOSAF	FWM2615
1.5cm	E7-1415 ISOSAF	FWM2622
1.7cm	E7-1417 ISOSAF	FWM2618
2.0cm	E7-1420 ISOSAF	FWM2616
2.3cm	E7-1423 ISOSAF	FWM2633
2.5cm	E7-1425 ISOSAF	FWM2613
2.7cm	E7-1427 ISOSAF	FWM2628
3.0cm	E7-1430 ISOSAF	FWM2617
3.5cm	E7-1435 ISOSAF	FWM2629
4.0cm	E7-1440 ISOSAF	FWM2623
4.4cm	E7-1444 ISOSAF	FWM2603
5.0cm	E7-1450 ISOSAF	FWM2570
5.5cm	E7-1455 ISOSAF	FWM2584
6.0cm	E7-1460 ISOSAF	FWM2579
6.5cm	E7-1465 ISOSAF	FWM2580

16Fr	Direct Code	NHS SC Code
1.0cm	E7-1610 ISOSAF	FWM2599
1.2cm	E7-1612 ISOSAF	FWM2589
1.5cm	E7-1615 ISOSAF	FWM2571
1.7cm	E7-1617 ISOSAF	FWM2619
2.0cm	E7-1620 ISOSAF	FWM2630
2.3cm	E7-1623 ISOSAF	FWM2585
2.5cm	E7-1625 ISOSAF	FWM2595
2.7cm	E7-1627 ISOSAF	FWM2572
3.0cm	E7-1630 ISOSAF	FWM2586
3.5cm	E7-1635 ISOSAF	FWM2581
4.0cm	E7-1640 ISOSAF	FWM2568
4.4cm	E7-1644 ISOSAF	FWM2590

18Fr	Direct Code	NHS SC Code
1.0cm	E7-1810 ISOSAF	FWM2689
1.2cm	E7-1812 ISOSAF	FWM2659
1.5cm	E7-1815 ISOSAF	FWM2569
1.7cm	E7-1817 ISOSAF	FWM2625
2.0cm	E7-1820 ISOSAF	FWM2624
2.3cm	E7-1823 ISOSAF	FWM2591
2.5cm	E7-1825 ISOSAF	FWM2627
2.7cm	E7-1827 ISOSAF	FWM2701
3.0cm	E7-1830 ISOSAF	FWM2597
3.5cm	E7-1835 ISOSAF	FWM2604
4.0cm	E7-1840 ISOSAF	FWM2601
4.4cm	E7-1844 ISOSAF	FWM2709

20Fr	Direct Code	NHS SC Code
1.0cm	E7-2010 ISOSAF	N/A
1.2cm	E7-2012 ISOSAF	FWM2710
1.5cm	E7-2015 ISOSAF	FWM2575
1.7cm	E7-2018 ISOSAF	FWM2667
2.0cm	E7-2020 ISOSAF	FWM2592
2.3cm	E7-2023 ISOSAF	N/A
2.5cm	E7-2025 ISOSAF	FWM2582
2.7cm	E7-2027 ISOSAF	FWM2702
3.0cm	E7-2030 ISOSAF	FWM2605
3.5cm	E7-2035 ISOSAF	FWM2573
4.0cm	E7-2040 ISOSAF	FWM2796
4.4cm	E7-2044 ISOSAF	FWM2711

24Fr	Direct Code	NHS SC Code
1.5cm	E7-2415 ISOSAF	FWM2574
1.7cm	E7-2417 ISOSAF	FWM2660
2.0cm	E7-2420 ISOSAF	FWM2587
2.3cm	E7-2423 ISOSAF	FWM2593
2.5cm	E7-2425 ISOSAF	N/A
2.7cm	E7-2427 ISOSAF	FWM2696
3.0cm	E7-2430 ISOSAF	FWM2703
3.5cm	E7-2435 ISOSAF	FWM2668
4.0cm	E7-2440 ISOSAF	N/A
4.4cm	E7-2444 ISOSAF	FWM2677

**Non-stock item:** orders placed on request available in 4-6 weeks.



### What is an E7?

The Pack includes: A single MINI Balloon Button E7 with a 5ml syringe G-Pack 1x12" right angle ENFit feeding set.



SM1-5  
ISO 5837

Mini  
ONE<sup>+</sup>  
BALLOON  
BUTTON

### What is an SM1-5?

A single MiniONE Balloon Button with 5ml syringe & introducer.

GLOW IN THE DARK PORT



## MiniONE Balloon Button SM1-5 Order Codes

12Fr	Direct Code	NHS SC Code
0.8cm	SM1-5-1208	FWM3046
1.0cm	SM1-5-1210	FWM3047
1.2cm	SM1-5-1212	FWM3048
1.5cm	SM1-5-1215	FWM3049
1.7cm	SM1-5-1217	FWM3050
2.0cm	SM1-5-1220	FWM3051
2.3cm	SM1-5-1223	FWM3052
2.5cm	SM1-5-1225	FWM3053
2.7cm	SM1-5-1227	FWM3054
3.0cm	SM1-5-1230	FWM3055
3.5cm	SM1-5-1235	FWM3056
4.0cm	SM1-5-1240	FWM3057

14Fr	Direct Code	NHS SC Code
0.8cm	SM1-5-1408	FWM3058
1.0cm	SM1-5-1410	FWM3059
1.2cm	SM1-5-1412	FWM3060
1.5cm	SM1-5-1415	FWM3061
1.7cm	SM1-5-1417	FWM3062
2.0cm	SM1-5-1420	FWM3063
2.3cm	SM1-5-1423	FWM3064
2.5cm	SM1-5-1425	FWM3065
2.7cm	SM1-5-1427	FWM3066
3.0cm	SM1-5-1430	FWM3067
3.5cm	SM1-5-1435	FWM3068
4.0cm	SM1-5-1440	FWM3069
4.4cm	SM1-5-1444	FWM3070
5.0cm	SM1-5-1450	FWM3071
5.5cm	SM1-5-1455	FWM3072
6.0cm	SM1-5-1460	FWM3073
6.5cm	SM1-5-1465	FWM3074

16Fr	Direct Code	NHS SC Code
1.0cm	SM1-5-1610	FWM3075
1.2cm	SM1-5-1612	FWM3076
1.5cm	SM1-5-1615	FWM3077
1.7cm	SM1-5-1617	FWM3078
2.0cm	SM1-5-1620	FWM3079
2.3cm	SM1-5-1623	FWM3080
2.5cm	SM1-5-1625	FWM3081
2.7cm	SM1-5-1627	FWM3082
3.0cm	SM1-5-1630	FWM3083
3.5cm	SM1-5-1635	FWM3084
4.0cm	SM1-5-1640	FWM3085
4.4cm	SM1-5-1644	FWM3086

18Fr	Direct Code	NHS SC Code
1.0cm	SM1-5-1810	FWM3087
1.2cm	SM1-5-1812	FWM3088
1.5cm	SM1-5-1815	FWM3089
1.7cm	SM1-5-1817	FWM3090
2.0cm	SM1-5-1820	FWM3091
2.3cm	SM1-5-1823	FWM3092
2.5cm	SM1-5-1825	FWM3093
2.7cm	SM1-5-1827	FWM3094
3.0cm	SM1-5-1830	FWM3095
3.5cm	SM1-5-1835	FWM3096
4.0cm	SM1-5-1840	FWM3097
4.4cm	SM1-5-1844	FWM3098

20Fr	Direct Code	NHS SC Code
1.0cm	SM1-5-2010	FWM3099
1.2cm	SM1-5-2012	FWM3100
1.5cm	SM1-5-2015	FWM3101
1.7cm	SM1-5-2018	FWM3102
2.0cm	SM1-5-2020	FWM3103
2.3cm	SM1-5-2023	FWM3104
2.5cm	SM1-5-2025	FWM3105
2.7cm	SM1-5-2027	FWM3106
3.0cm	SM1-5-2030	FWM3107
3.5cm	SM1-5-2035	FWM3108
4.0cm	SM1-5-2040	FWM3109
4.4cm	SM1-5-2044	FWM3110

24Fr	Direct Code	NHS SC Code
1.5cm	SM1-5-2415	FWM3111
1.7cm	SM1-5-2417	FWM3112
2.0cm	SM1-5-2420	FWM3113
2.3cm	SM1-5-2423	FWM3114
2.5cm	SM1-5-2425	FWM3115
2.7cm	SM1-5-2427	FWM3116
3.0cm	SM1-5-2430	FWM3117
3.5cm	SM1-5-2435	FWM3118
4.0cm	SM1-5-2440	FWM3119
4.4cm	SM1-5-2444	FWM3120

MiniONE Balloon Button Feeding Sets

See page 25

Non-stock item: orders placed on request available in 4-6 weeks.

## MiniONE Balloon Button E7-M1 Order Codes

12Fr	Direct Code	NHS SC Code
0.8cm	E7-M1-1208 ISOSAF	N/A
1.0cm	E7-M1-1210 ISOSAF	N/A
1.2cm	E7-M1-1212 ISOSAF	N/A
1.5cm	E7-M1-1215 ISOSAF	FWM3128
1.7cm	E7-M1-1217 ISOSAF	FWM3129
2.0cm	E7-M1-1220 ISOSAF	FWM3130
2.3cm	E7-M1-1223 ISOSAF	FWM3131
2.5cm	E7-M1-1225 ISOSAF	FWM3132
2.7cm	E7-M1-1227 ISOSAF	FWM3133
3.0cm	E7-M1-1230 ISOSAF	FWM3134
3.5cm	E7-M1-1235 ISOSAF	FWM3135
4.0cm	E7-M1-1240 ISOSAF	FWM3136

14Fr	Direct Code	NHS SC Code
0.8cm	E7-M1-1408 ISOSAF	FWM3137
1.0cm	E7-M1-1410 ISOSAF	FWM3138
1.2cm	E7-M1-1412 ISOSAF	FWM3139
1.5cm	E7-M1-1415 ISOSAF	FWM3140
1.7cm	E7-M1-1417 ISOSAF	FWM3141
2.0cm	E7-M1-1420 ISOSAF	FWM3142
2.3cm	E7-M1-1423 ISOSAF	FWM3143
2.5cm	E7-M1-1425 ISOSAF	FWM3144
2.7cm	E7-M1-1427 ISOSAF	FWM3145
3.0cm	E7-M1-1430 ISOSAF	FWM3146
3.5cm	E7-M1-1435 ISOSAF	FWM3147
4.0cm	E7-M1-1440 ISOSAF	FWM3148
4.4cm	E7-M1-1444 ISOSAF	N/A
5.0cm	E7-M1-1450 ISOSAF	N/A
5.5cm	E7-M1-1455 ISOSAF	FWM3149
6.0cm	E7-M1-1460 ISOSAF	N/A
6.5cm	E7-M1-1465 ISOSAF	N/A

16Fr	Direct Code	NHS SC Code
1.0cm	E7-M1-1610 ISOSAF	N/A
1.2cm	E7-M1-1612 ISOSAF	N/A
1.5cm	E7-M1-1615 ISOSAF	N/A
1.7cm	E7-M1-1617 ISOSAF	FWM3150
2.0cm	E7-M1-1620 ISOSAF	N/A
2.3cm	E7-M1-1623 ISOSAF	N/A
2.5cm	E7-M1-1625 ISOSAF	N/A
2.7cm	E7-M1-1627 ISOSAF	N/A
3.0cm	E7-M1-1630 ISOSAF	N/A
3.5cm	E7-M1-1635 ISOSAF	N/A
4.0cm	E7-M1-1640 ISOSAF	N/A
4.4cm	E7-M1-1644 ISOSAF	N/A

18Fr	Direct Code	NHS SC Code
1.0cm	E7-M1-1810 ISOSAF	N/A
1.2cm	E7-M1-1812 ISOSAF	N/A
1.7cm	E7-M1-1817 ISOSAF	FWM3151
2.0cm	E7-M1-1820 ISOSAF	FWM3152
2.3cm	E7-M1-1823 ISOSAF	FWM3153
2.5cm	E7-M1-1825 ISOSAF	FWM3154
2.7cm	E7-M1-1827 ISOSAF	N/A
3.0cm	E7-M1-1830 ISOSAF	FWM3155
3.5cm	E7-M1-1835 ISOSAF	FWM3156
4.0cm	E7-M1-1840 ISOSAF	N/A
4.4cm	E7-M1-1844 ISOSAF	FWM3157

20Fr	Direct Code	NHS SC Code
1.0cm	E7-M1-2010 ISOSAF	N/A
1.2cm	E7-M1-2012 ISOSAF	N/A
1.7cm	E7-M1-2018 ISOSAF	N/A
2.0cm	E7-M1-2020 ISOSAF	FWM3158
2.3cm	E7-M1-2023 ISOSAF	NA
2.5cm	E7-M1-2025 ISOSAF	FWM3159
2.7cm	E7-M1-2027 ISOSAF	NA
3.0cm	E7-M1-2030 ISOSAF	FWM3160
3.5cm	E7-M1-2035 ISOSAF	FWM3161
4.0cm	E7-M1-2040 ISOSAF	FWM3162

24Fr	Direct Code	NHS SC Code
1.7cm	E7-M1-2417 ISOSAF	N/A
2.0cm	E7-M1-2420 ISOSAF	N/A
2.3cm	E7-M1-2423 ISOSAF	N/A
2.5cm	E7-M1-2425 ISOSAF	N/A
2.7cm	E7-M1-2427 ISOSAF	N/A
3.0cm	E7-M1-2430 ISOSAF	N/A
3.5cm	E7-M1-2435 ISOSAF	N/A
4.0cm	E7-M1-2440 ISOSAF	N/A

	Direct Code	Product Description	NHS SC Code
12Fr	E7-M1-1208 ISOSAF	Mini One balloon button 12Fr 0.8cm with syringe and introducer Gpak and 1 x 12" right and ENFit feeding set	FWM3261
12Fr	E7-M1-1210 ISOSAF	Mini One balloon button 12Fr 1.0cm with syringe and introducer Gpak and 1 x 12" right and ENFit feeding set	FWM3262
12Fr	E7-M1-1212 ISOSAF	Mini One balloon button 12Fr 1.2cm with syringe and introducer Gpak and 1 x 12" right and ENFit feeding set	FWM3263
14Fr	E7-M1-1444 ISOSAF	Mini One balloon button 14Fr 4.4cm with syringe and introducer Gpak and 1 x 12" right and ENFit feeding set	FWM3264
14Fr	E7-M1-1450 ISOSAF	Mini One balloon button 14Fr 5.0cm with syringe and introducer Gpak and 1 x 12" right and ENFit feeding set	FWM3265
14Fr	E7-M1-1460 ISOSAF	Mini One balloon button 14Fr 6.0cm with syringe and introducer Gpak and 1 x 12" right and ENFit feeding set	FWM3266



### What is an E7-M1?

The Pack includes:  
A single MiniONE Balloon Button with a 5ml syringe G-Pak 1x 12" right angle ENFit feeding set.



8-1223-ISOSAF

**Non-stock item:** orders placed on request available in 4-6 weeks.

**G-PAK**  
ISO Saf

### G-Pak

The pack contains everything you need for a Gastrostomy change in one easy-to-use pack.



Shown are for illustration purposes only. Actual product may vary.

### G-Pak Order Codes

Direct Code	NHS SC Code	Unit of Sale
G-PAKISO	FWM2654	1

### Kit Contains

QTY	Wrapper
1	Stericrepe Paper 60 x 60cm
<b>Layer 1</b>	
2	Gloves - Latex Free
1	Plastic Apron
1	Lubrication Aqua Gel
1x5	Swab, Plain, Gauze 10 x 10cm, 12-ply White
1	Syringe 10ml Sterile Water Prefilled
1	Syringe, 10ml Luer Slip Eccentric
1	Syringe, 60ml ENFit ISO 80369-3
1	Tray Gallipot
	Waste Bag

### Which Accessories & Related Products Can I Use?

**BALLOON G-TUBE**  
ISO Saf



See page 16

**CAPSULE MONARCH NON-BALLOON G-TUBE**  
ISO Saf



See page 17

**MINI BALLOON BUTTON**



See page 18-19

**G-JET LOW PROFILE GASTRO-JEJUNAL ENTERAL TUBE**



See page 26-27

## MiniONE Non-Balloon Button

MiniONE Non-Balloon Button lies flush against the skin of the abdomen and has a small internal bolster that sits inside the stomach to prevent the device from falling out.



### MiniONE Non-Balloon Button Order Codes

14Fr	Direct Code	NHS SC Code
1.0cm	M1-2-1410-I	FWM2598
1.2cm	M1-2-1412-I	FWM1393
1.5cm	M1-2-1415-I	FWM2602
1.7cm	M1-2-1417-I	FWM2600
2.0cm	M1-2-1420-I	FWM2576
2.5cm	M1-2-1425-I	FWM2594
3.0cm	M1-2-1430-I	FWM2583

18Fr	Direct Code	NHS SC Code
3.4cm	M1-2-1834-I	FWM2669
4.4cm	M1-2-1844-I	FWM2712

### MiniONE Non-Balloon Button Removal Devices

14Fr	Direct Code	NHS SC Code
1.0-1.2cm	R-S1012	N/A
1.5-1.7cm	R-S1517	FWM1362
2.0cm	R-S2023	FWM1359
2.5cm	R-S2427	FWM1360
3.0cm	R-S3033	N/A

18Fr	Direct Code	NHS SC Code
3.4-3.7cm	R-L3437	N/A
4.0-4.4cm	R-L4044	N/A

### Shaft Length (cm)

	1.0	1.2	1.5	1.7	2.0	2.5	3.0	3.4	4.4
14Fr	●	●	●	●	●	●	●		
18Fr								●	●

MINI  
BALLOON BUTTON













FEED  
SETS  
ISOsaf

Sets can be used with...

MINI  
BALLOON BUTTON

G-JET<sup>†</sup>  
LOW PROFILE  
GASTRO-JEJUNAL  
ENTERAL TUBE

MINI Balloon Button Feed Set Order Codes

	Size	Description	Direct Code	NHS SC Code	Unit of Sale
	2"	ENFit Y-Port, Straight	6-0212 ISOSAF	FWM2796	10
	12"	ENFit Port, Right angle	6-1223 ISOSAF	FWM2425	10
	24"	ENFit Port, Right angle	6-2423 ISOSAF	FWM2412	10
	12"	ENFit Port, Straight with tethered cap	6-1211 ISOSAF	FWM2424	10
	24"	ENFit Port, Straight with tethered cap	6-2411 ISOSAF	FWM2723	10
	12"	ENFit Port, Right angle with tethered cap	6-1226 ISOSAF	FWM2411	10
	24"	ENFit Port, Right angle with tethered cap	6-2421 ISOSAF	N/A	10
	12"	ENFit Port, Right angle with tethered cap, sterile	S6-1226 ISOSAF	FWM2642	10
	12"	ENFit Port, ENFit med port with right angle	6-1224 ISOSAF	FWM2722	10
	24"	ENFit Port, ENFit med port with right angle	6-2424 ISOSAF	FWM2721	10
	12"	ENFit Port, Right Angle, Double lock, Single	6-1223D ISOSAF	FWM2636	10
	24"	ENFit Port, Right Angle, Double lock, Single	6-2423D ISOSAF	FWM2606	10












6-Feeding Sets are also used for Gastrostomy venting with the **AMT G-Jet**




Sets can be used with...



### MiniONE Feed Set Order Codes

	Size	Description	Direct Code	NHS SC Code	Unit of Sale
	2"	ENFit Y-Port, Straight	8-0212 ISOSAF	FWM2717	10
	12"	ENFit Port, Right angle	8-1223 ISOSAF	FWM2715	10
	24"	ENFit Port, Right angle	8-2423 ISOSAF	FWM2609	10
	12"	ENFit Port, Straight with tethered cap	8-1211 ISOSAF	FWM2714	10
	24"	ENFit Port, Straight with tethered cap	8-2411 ISOSAF	FWM2720	10
	12"	ENFit Port, Right angle with tethered cap	8-1226 ISOSAF	FWM2611	10
	12"	ENFit Port, ENFit med port with right angle	8-1224 ISOSAF	FWM2719	10
	24"	ENFit Port, ENFit med port with right angle	8-2424 ISOSAF	FWM2718	10
	12"	ENFit Port, Right angle double lock	8-1223D ISOSAF	FWM2713	10
	24"	ENFit Port, Right angle double lock	8-2423D ISOSAF	FWM2716	10
	12"	ENFit Y Port, Glow in the dark right angle	8-1255 ISOSAF	FWM2612	10

8-Feeding Sets are also used for Jejunal feeding with the **AMT G-Jet**



**Retention Balloon**

- In order to feed, give medication, or vent from the AMT G-Jet, a feeding/venting set needs to be attached.

**AMT G-Jet**

The AMT G-Jet is a triple lumen feeding device which is used to help patients meet their nutritional needs by feeding directly into the small bowel whilst being able to vent the stomach at the same time.

G-Jet Order Codes

14Fr		Direct Code	Availability
Gastric	Jejunal		
1.0cm	15cm	GJ-1410-15-I	Stock Item
1.0cm	22cm	GJ-1410-22-I	Stock Item
1.2cm	15cm	GJ-1412-15-I	Stock Item
1.2cm	22cm	GJ-1412-22-I	Stock Item
1.2cm	30cm	GJ-1412-30-I	Stock Item
1.5cm	15cm	GJ-1415-15-I	4-6 Weeks
1.5cm	22cm	GJ-1415-22-I	Stock Item
1.5cm	30cm	GJ-1415-30-I	Stock Item
1.5cm	45cm	GJ-1415-45-I	4-6 Weeks
1.7cm	15cm	GJ-1417-15-I	4-6 Weeks
1.7cm	22cm	GJ-1417-22-I	Stock Item
1.7cm	30cm	GJ-1417-30-I	Stock Item
1.7cm	45cm	GJ-1417-45-I	4-6 Weeks
2.0cm	22cm	GJ-1420-22-I	Stock Item
2.0cm	30cm	GJ-1420-30-I	Stock Item
2.0cm	45cm	GJ-1420-45-I	Stock Item

14Fr		Direct Code	Availability
Gastric	Jejunal		
2.3cm	30cm	GJ-1423-30-I	Stock Item
2.3cm	45cm	GJ-1423-45-I	Stock Item
2.5cm	30cm	GJ-1425-30-I	Stock Item
2.5cm	45cm	GJ-1425-45-I	Stock Item
2.7cm	30cm	GJ-1427-30-I	Stock Item
2.7cm	45cm	GJ-1427-45-I	4-6 Weeks
3.0cm	45cm	GJ-1430-45-I	Stock Item
3.5cm	45cm	GJ-1435-45-I	Stock Item
4.0cm	45cm	GJ-1440-45-I	Stock Item
4.4cm	45cm	GJ-1444-45-I	Stock Item
5.0cm	45cm	GJ-1450-45-I	Stock Item
5.5cm	45cm	GJ-1455-45-I	Stock Item
6.0cm	45cm	GJ-1460-45-I	Stock Item

**Non-stock item:** orders placed on request available in 4-6 weeks.



**Gastric Venting Sets** See page 24

- The code on the packaging has to start with **6-xxxx**.
- The end which attaches to the G-port should be white.

**Balloon Inflation Port**

- There are a variety available to suit patients' needs.
- If both ports are used at the same time, attach the venting set into the G-port first.

**Jejunal Feeding Sets** See page 25

- The code on the packaging has to start with **8-xxxx**.
- The end which attaches to the J-port should be beige or luminous.

16Fr		Direct Code	Availability
Gastric	Jejunal		
1.0cm	15cm	GJ-1610-15-I	4-6 Weeks
1.0cm	22cm	GJ-1610-22-I	4-6 Weeks
1.2cm	15cm	GJ-1612-15-I	4-6 Weeks
1.2cm	22cm	GJ-1612-22-I	Stock Item
1.2cm	30cm	GJ-1612-30-I	4-6 Weeks
1.5cm	15cm	GJ-1615-15-I	4-6 Weeks
1.5cm	22cm	GJ-1615-22-I	4-6 Weeks
1.5cm	30cm	GJ-1615-30-I	Stock Item
1.5cm	45cm	GJ-1615-45-I	4-6 Weeks
1.7cm	15cm	GJ-1617-15-I	Stock Item
1.7cm	22cm	GJ-1617-22-I	4-6 Weeks
1.7cm	30cm	GJ-1617-30-I	Stock Item
1.7cm	45cm	GJ-1617-45-I	Stock Item
2.0cm	22cm	GJ-1620-22-I	Stock Item
2.0cm	30cm	GJ-1620-30-I	Stock Item

16Fr		Direct Code	Availability
Gastric	Jejunal		
2.0cm	45cm	GJ-1620-45-I	Stock Item
2.3cm	30cm	GJ-1623-30-I	Stock Item
2.3cm	45cm	GJ-1623-45-I	Stock Item
2.5cm	30cm	GJ-1625-30-I	Stock Item
2.5cm	45cm	GJ-1625-45-I	Stock Item
2.7cm	30cm	GJ-1627-30-I	Stock Item
2.7cm	45cm	GJ-1627-45-I	Stock Item
3.0cm	45cm	GJ-1630-45-I	Stock Item
3.5cm	45cm	GJ-1635-45-I	Stock Item


**Non-stock item:** orders placed on request available in 4-6 weeks.

**ENPLUG**

**ENPLUG Stoma Plugs**

ENPLUG are designed for emergency use to keep stomas from closing when a gastrostomy device has fallen out.

The ENPLUG packs are available in 2 different lengths and various French sizes in case the stoma has already started to close and the current French gauge tube no longer fits the stoma.



ENPLUG Order Codes

Size	Pack Content	Codes	NHS SC Code	Length	Packaging
S	10,12,14,16Fr	EN-PLUG-S	FWM3197	4cm	4pcs/pouch
L	10,12,14,16Fr	EN-PLUG-L	FWM3198	7cm	4pcs/pouch
XL	18, 20Fr	EN-PLUG-XL	FWM3199	7cm	2pcs/pouch

Which Accessories & Related Products Can I Use?

**CAPSULE MONARCH**  
NON-BALLOON  
G-TUBE  
ISO 5834

See page 17

**MINI**  
BALLOON  
BUTTON

See page 18-19

**MINI ONE**  
NON-BALLOON  
BUTTON

See page 23

**G-JET**  
LOW PROFILE  
GASTRO-JEJUNAL  
ENTERAL TUBE

See page 26-27



**ENTERAL  
SYRINGE  
DRIVER**

## Enteral Syringe Driver Order Codes

Direct Code	NHS SC Code
ENT575	N/A

## Which Accessories & Related Products Can I Use?



**ENTERAL**  
ISO SAF

## Enteral Single Use Syringes

Enteral ISO SAF, single use enteral syringes with female ENFit connection (ISO 80369-3). Low dose option available in 1ml & 3ml syringes.



### Enteral ISO SAF Order Codes

Size	Low Dose	Codes	NHS SC Code	Packaging	Unit of Sale
1ml	•	ELDISO-1	FWM2427	1pc/pouch	100pcs./box
3ml	•	ELDISO-3	FWM2428	1pc/pouch	100pcs./box
1ml		EISO-1	FTA350	1pc/pouch	100pcs./box
3ml		EISO-3	FTA351	1pc/pouch	100pcs./box
5ml		EISO-5	FWM2430	1pc/pouch	100pcs./box
10ml		EISO-10	FWM2431	1pc/pouch	100pcs./box
20ml		EISO-20	FWM2426	1pc/pouch	100pcs./box
60ml		EISO-60	FWM2429	1pc/pouch	50pcs./box

### Product Features



### Which Accessories & Related Products Can I Use?



## Enteral Home Syringes

Enteral ISO SAF Home, reusable enteral syringe with female ENFit connection (ISO 80369-3). Low dose connections are fitted to the 1ml and 3ml syringes.



**ENTERAL**  
ISO SAF Home

LOW DOSE  
SYRINGE

### Enteral ISO SAF Home Order Codes

Size	Low Dose	Codes	NHS SC Code	Packaging	Unit of Sale
1ml	•	ELDISO-1H	FTR2105	1pc/pouch	100pcs./box
3ml	•	ELDISO-3H	FTR2102	1pc/pouch	100pcs./box
5ml		EISO-5H	FTR2101	1pc/pouch	100pcs./box
10ml		EISO-10H	FTR2103	1pc/pouch	100pcs./box
20ml		EISO-20H	FTR2104	1pc/pouch	100pcs./box
60ml		EISO-60H	FTR2106	1pc/pouch	50pcs./box

### Product Features

FULL  
THREAD

LATEX  
FREE

DEHP  
FREE

UP TO  
**84**  
USES

RUBBER  
GASKET

### Which Accessories & Related Products Can I Use?



CareTip  
**ORAL**  
SYRINGE




Oral Syringe Order Codes

Size	Codes	NHS SC Code	Packaging	Unit of Sale
1ml	OS-1	FTA343	1pc/pouch	100pcs./box
2.5ml	OS-2	FTA344	1pc/pouch	100pcs./box
5ml	OS-5	FTA345	1pc/pouch	100pcs./box
10ml	OS-10	FTA355	1pc/pouch	100pcs./box

Product Features



How To Use An Oral Syringe



Related Products







**CareTip**  
**CATHETER TIP**  
**SYRINGE**

### Catheter Tip Order Code

Length	Codes	NHS SC Code	Packaging	Unit of Sale
60ml	CT60	FWC449	1pc./pouch	50pcs./box

### Product Features

**DEHP**  
**FREE**

**LATEX**  
**FREE**

### Colostrum Gold Kit

Kit for the collection and administration of Colostrum available with 1ml or 3ml syringes.

#### Kit Contains

QTY	Product Description
1	Information Leaflet
2	Bonding Heart
10	Identification Label
10	Syringe, 1ml or 3ml
10	Syringe Cap



**COLOSTRUM**  
**GOLD KIT**

### Colostrum Gold Kit Order Codes

Size	Codes	NHS SC Code	Packaging	Unit of Sale
1ml	CGK-1	FTA366	1pc./pouch	10pcs./box
3ml	CGK-3	FTA367	1pc./pouch	10pcs./box

ISO Saf  
ACCESSORIES

ENFit  
BOTTLE  
ADAPTOR  
ISO Saf

ENFit Bottle Adaptor

Available in five sizes to fit the majority of medicine bottles.



Order Codes - Singles

Size	Direct Codes	NHS SC Code	Unit Of Sale
12 – 15mm	BAISO1-S	FDC1366	Single blister 100 / box
14 – 16.5mm	BAISO2-S	FDC1367	Single blister 100 / box
17– 20.5mm	BAISO3-S	FDC1368	Single blister 100 / box
19 - 23mm	BAISO35-S	FDC1369	Single blister 100 / box
24 – 26mm	BAISO4-S	FDC1370	Single blister 100 / box

Order Codes - Bulk Pack

Size	Direct Codes	NHS SC Code	Unit Of Sale
12 – 15mm	BAISO1	FDC1135	100pcs.
14 – 16.5mm	BAISO2	FDC1136	100pcs.
17– 20.5mm	BAISO3	FDC1137	100pcs.
19 - 23mm	BAISO35	FDC1349	100pcs.
24 – 26mm	BAISO4	FDC1138	100pcs.

ENTERAL  
MEDICINE  
STRAW  
ISO Saf

Enteral Medicine Straw

Male ENFit connector (ISO 80369-3).



Order Codes

Non Sterile

Size	Direct Code	NHS SC Code	Packaging	Unit Of Sale
50mm	MSISO-50	FWM2610	1pc/pouch	100 per box

Sterile

50mm	13051050	FDC1146	1pc/pouch	30 per box
100mm	13051100	FDC1123	1pc/pouch	30 per box
150mm	13051150	FDC1145	1pc/pouch	30 per box

MEDICINE  
BOTTLE  
ADAPTOR  
ISO Saf

Medicine Bottle Adaptor

Stepped Bottle Adaptor, 15-27mm, male ENFit connector (ISO 80369-3).



Order Codes

Direct Code	NHS SC Code	Packaging	Unit Of Sale
SBAISO	FDC1134	1pc/pouch	50 per box

ENFit®  
ACCESSORIES



### Enteral Syringe Cap

Male ENFit connector (ISO 80369-3), sterile.



#### Order Codes

Direct Code	NHS SC Code	Packaging	Unit Of Sale
ECISO	FWD276	1pc/pouch	100 per box
ECISO-10	FWD284	10pc/pouch	100 per box



### Filter Needle

Male ENFit connector (ISO 80369-3), sterile. 5µm.



#### Order Codes

Dash3				
Size	Direct Code	NHS SC Code	Packaging	Unit Of Sale
18G 38mm 5µm	13061038	FDC1147	1pc/pouch	30pcs./box
18G 50mm 5µm	13061050	FDC1148	1pc/pouch	30pcs./box



### ENFit Female Cap

Female ENFit connector (ISO 80369-3).



#### Order Codes

Direct Code	NHS SC Code	Packaging	Unit Of Sale
FCISO-NS	FWD277	200pc/pouch	200pcs./box
FCISO-10	FWD283	10pc/pouch	150pcs./box



### Y-Port Adaptor

Double male ENFit connector (ISO 80369-3), converts any ENFit single port feeding tube to Y-Port, sterile.



#### Order Codes

Direct Code	NHS SC Code	Packaging	Unit Of Sale
YPA-ISO	FWM2432	1pc/pouch	20 per box

ENTERAL  
ACCESSORIES

BUTTON  
BEAR



Order Codes

Colour	Direct Code
Blonde	P-Bear(B)
Brown	P-Bear(C)

FIXATION  
PLASTER

Fixation Plaster



Order Codes

Direct Code	NHS SC Code	Unit Of Sale
PLNC	N/A	100 pcs./box

ORAL  
BOTTLE  
ADAPTOR

Oral Bottle  
Adaptor



Order Codes

Size	Direct Code	NHS SC Code	Unit Of Sale
12mm-16mm	OBA-A	FDC1350	100pcs.
17mm-20mm	OBA-B	FDC1351	100pcs.
21mm-23mm	OBA-C	FDC1352	100pcs.
24mm-28mm	OBA-D	FDC1353	100pcs.

ENTERAL  
Halo

Enteral Halo  
Skin Protector



Order Codes

	Direct Code	NHS SC Code	Unit Of Sale
58m	ENT-HALO	FWM1890	10 pcs./box

ENTERAL  
CARE  
CASE



Order Codes

Direct Code
ECK002

ENTERAL  
DEMO  
CASE



Order Codes

Direct Code
ECK001


SIDE  
BOARD

Side Board

Training Board for NG  
Tube & Gastrostomy  
placement.

**Includes:**

- 1x NG tube
- 1x ENFit syringe
- 1x Silicone oil
- 1x Shoulder bag
- 2x Side brackets



Side Board Order Codes

Direct Code	NHS SC Code
TRSB01	N/A

ENFit® TO  
FUNNEL  
ADAPTOR



Order Codes

Direct Code	NHS SC Code	Unit Of Sale
UAD2ISO	FWM2437	100 pcs./box

FUNNEL  
TO ENFit®  
ADAPTOR



Order Codes

Direct Code	NHS SC Code	Unit Of Sale
TRN202	FWM2436	10 pcs./box

AMT  
ACCESSORIES

AMT  
INTRODUCER

AMT Introducer



Order Codes

Direct Code	NHS SC Code	Unit Of Sale
4-5000	FWM864	10 pcs./box

AMT  
STIFFENER  
GUIDEWIRE  
COMPATIBLE†

AMT Stiffener Guidewire  
Compatible



Order Codes

Direct Code	NHS SC Code	Unit Of Sale
4-6000	FWM2153	10 pcs./box

AMT  
STOMA  
MEASURING  
DEVICE†

AMT Stoma Measuring  
Device



Order Codes

Direct Code	NHS SC Code	Unit Of Sale
1704	FWM644	10 pcs./box

AMT  
STOMA  
MEASURING  
DEVICE  
GUIDEWIRE  
COMPATIBLE†

AMT Stoma Measuring Device  
Guidewire Compatible



Order Codes

Direct Code	NHS SC Code	Unit Of Sale
1804	FWM2152	10 pcs./box

AMT  
CLAMP†

AMT Clamp



Order Codes

Direct Code	NHS SC Code	Unit Of Sale
4-3000	FMW708	10 pcs./box




# Quality Everyday Essentials

Scan the QR codes to view brochures




## Critical Care

**Open Suction, Closed Suction,  
Blood Pressure Cuffs and CareTip Yankauers**




## Banana Patient Moving & Handling Brochure

**Air Assisted Devices  
Care Management Systems & Accessories  
Patient Repositioning, Proning and Transfer Devices**




## Blood Pressure Cuff Brochure

**Reusable BP Cuffs  
Single Patient Use BP Cuffs  
and Neonatal BP Cuffs**




## Serres Brochure

**Serres Suction Liner System**



## Wound Drainage Brochure

**High Vacuum Systems and Bellows Systems  
Silicone Systems and Robinson Type Systems**



Scan the QR  
code using your  
smart phone camera

  
**homesafe**  
direct to your door

  
**homesafe**  
direct to your door

**HEALTHCARE  
PROFESSIONAL  
BENEFITS**

- Most cost effective way to get Enteral UK products to your patients
- Fast, easy set up process
- Standard next working day delivery
- Secure data exchange through NHS.net
- Control of ordering: standing order, call off, ad hoc or a mix you decide
- Your invoices, your way
- Dedicated customer service contact
- Local Enteral Specialist support

**PATIENT  
BENEFITS  
FROM**

- Tailored welcome package
- Customer helpline
- Call off & emergency orders delivered next working day
- Deliveries synchronised with enteral feeding company
- Minimal packaging, reduced space needed for storage
- Exchange programme for size changes
- Multilingual product information



 **Enteral**  
UK

  
Applied Medical Technology Inc.

GBUK Enteral Ltd. Woodland House, Blackwood Hall Business Park  
North Duffield, Selby, North Yorkshire, YO8 5DD, UK.

T: +44 (0)1757 282 945  
F: +44 (0)1757 600 545

info@gbukenteral.com  
www.gbukenteral.com

Registered Trademark or Trademark of GBUK Enteral Ltd. † Registered Trademark or Trademark of Applied Medical Technology Inc.  
ENFit® is a Trademark of GEDSA. MIC-KEY is a registered Trademark of Halyard Health, Inc.,

Enteral UK is a trading name of GBUK Enteral Ltd. ©2019 Copyright GBUK Group Ltd.

ENT PL002 

05/19