

# Export





Enteral UK is dedicated to the provision of innovative devices that enable healthcare professionals to provide enteral feeding support to their patients.

Enteral UK is working with the world's leading enteral feeding device suppliers to implement the new ENFit® ISO 80369-3 international standard. We are collaborating with clinicians and representatives from all clinical interest groups and regulatory agencies to ensure that the transition to the new connectors is easy and safe for the patient.



GBUK Group award-winning specialist companies are expert in the provision of enteral feeding, critical care and patient moving and handling products and technologies. GBUK Group companies are suppliers to every UK NHS Hospital and to a rapidly expanding Export customer base.

Since GBUK was formed in July 2008 the organization has continually delivered innovative products based on 'needs' identified from Clinicians' feedback regarding the health, safety and the well-being of their patients. We believe this is what makes the UK health industry the best in the world; UK Clinicians are prepared to get involved and give their precious time to help companies like ours design, innovate and manufacture leading-edge products. However, we are not working on expensive medical solutions, we are working at delivering innovation and change to routinely used 'day-in day-out' medical products. We believe that our heuristic approach is unique within the healthcare industry and has resulted in the usage of GBUK products rising from 10,000 per day in 2008 to more than 213,000 today.

We have an active research and development team permanently engaged in discovering and developing new innovative products, as well as improving our existing ones. This enables us to continually deliver effective solutions and ensure that GBUK remains at the forefront of healthcare innovation.

GBUK Enteral Export



### GBUK Enteral


is the leading specialist enteral feeding device company in the UK. Our innovative products enable healthcare professionals to provide enteral feeding support to their patients. Every day thousands of patients in hospital and their own homes benefit from our products.



# Enteral Syringes

## Single Use Syringe

Enteral ISOSAF Single Use enteral syringe with female ENFit connection (ISO 80369-3). Low dose connections are available.



### Product Features

- Low dose 1ml, 3ml & 6ml syringes
- 100ml syringe now available
- CE marked for Enteral/Oral
- CE for measurement
- Full thread
- Rubber gasket in single use
- Latex & DEHP free
- Clear barrel markings

ENFit®

Direct Codes	Specification	UOS
<b>ELDISO-PFC-1</b>	ISOSAF™ Enteral Syringe, Low Dose, Single Use, ISO 80369-3, 1ml	100 units
<b>ELDISO-PFC-3</b>	ISOSAF™ Enteral Syringe, Low Dose, Single Use, ISO 80369-3, 3ml	100 units
<b>ELDISO-PFC-6</b>	ISOSAF™ Enteral Syringe, Low Dose, Single Use, ISO 80369-3, 6ml	100 units
<b>EISO-PFC-12</b>	ISOSAF™ Enteral Syringe, Single Use, ISO 80369-3, 12ml	100 units
<b>EISO-PFC-20</b>	ISOSAF™ Enteral Syringe, Single Use, ISO 80369-3, 20ml	100 units
<b>EISO-PFC-35</b>	ISOSAF™ Enteral Syringe, Single Use, ISO 80369-3, 35ml	50 units
<b>EISO-PFC-60</b>	ISOSAF™ Enteral Syringe, Single Use, ISO 80369-3, 60ml	50 units
<b>EISO-100</b>	ISOSAF™ Enteral Syringe, Single Use, ISO 80369-3, 100ml	30 units
<b>SELDISO-1</b>	ISOSAF™ Enteral Syringe, Low Dose, Single O Ring Gasket, Single Use, ISO 80369-3, 1ml	100 units
<b>SELDISO-3</b>	ISOSAF™ Enteral Syringe, Low Dose, Single O Ring Gasket, Single Use, ISO 80369-3, 3ml	100 units
<b>SELDISO-6</b>	ISOSAF™ Enteral Syringe, Low Dose, Single O Ring Gasket, Single Use, ISO 80369-3, 6ml	100 units
<b>SEISO-12</b>	ISOSAF™ Enteral Syringe, Single O Ring Gasket, Single Use, ISO 80369-3, 12ml	100 units
<b>SEISO-20</b>	ISOSAF™ Enteral Syringe, Single O Ring Gasket, Single Use, ISO 80369-3, 20ml	100 units
<b>SEISO-35</b>	ISOSAF™ Enteral Syringe, Single O Ring Gasket, Single Use, ISO 80369-3, 35ml	100 units
<b>SEISO-60</b>	ISOSAF™ Enteral Syringe, Single O Ring Gasket, Single Use, ISO 80369-3, 60ml	50 units

# Enteral Syringes

## "O" Ring Gasket Single Use Syringe

Single use "O" ring design.  
Syringes available in sizes 1-60ml.



ENFit connection  
(ISO 80369-3).  
Low dose in  
1ml, 3ml & 6ml  
syringes.

### With Push Fit Tip Caps



### Product Features

- Low dose 1ml, 3ml & 6ml syringes
- "O" ring gasket
- Silicone gasket
- Full thread
- Latex & DEHP free
- Clear barrel markings
- 60ml syringe has an eccentric tip for optimization of fat delivery when delivering milk

ENFit®

Direct Codes	Product Descriptions	Packaging	Unit of Sale
<b>SELDISO-PFC-1</b>	ISOSAF™ Enteral Syringe, Low Dose, Single O Ring Gasket and Push Fit Cap, Single Use, ISO 80369-3, 1ml	1pc/ Pouch	100 units/box
<b>SELDISO-PFC-3</b>	ISOSAF™ Enteral Syringe, Low Dose, Single O Ring Gasket and Push Fit Cap, Single Use, ISO 80369-3, 3ml	1pc/ Pouch	100 units/box
<b>SELDISO-PFC-6</b>	ISOSAF™ Enteral Syringe, Low Dose, Single O Ring Gasket and Push Fit Cap, Single Use, ISO 80369-3, 6ml	1pc/ Pouch	100 units/box
<b>SEISO-PFC-12</b>	ISOSAF™ Enteral Syringe, Single O Ring Gasket and Push Fit Cap, Single Use, ISO 80369-3, 12ml	1pc/ Pouch	100 units/box
<b>SEISO-PFC-20</b>	ISOSAF™ Enteral Syringe, Single O Ring Gasket and Push Fit Cap, Single Use, ISO 80369-3, 20ml	1pc/ Pouch	100 units/box
<b>SEISO-PFC-35</b>	ISOSAF™ Enteral Syringe, Single O Ring Gasket and Push Fit Cap, Single Use, ISO 80369-3, 35ml	1pc/ Pouch	50 units/box
<b>SEISO-PFC-60</b>	ISOSAF™ Enteral Syringe, Single O Ring Gasket and Push Fit Cap, Single Use, ISO 80369-3, 60ml	1pc/ Pouch	50 units/box



# Amber ENFit® Syringes

## Amber Barrel Enteral ISOSAF Single Use Syringe

Enteral ISOSAF Single Use enteral syringe with female ENFit connection (ISO 80369-3).

Amber barrel syringes do not transmit more than 10% of incident light with wavelengths of 290-450 nm. As per USP671.



### Product Features


- Low dose 1ml & 3ml syringes
- Full thread
- Rubber gasket in single use
- Latex & DEHP free
- Clear barrel markings

Direct Codes	Product Descriptions	Packaging	Unit of Sale
<b>ABLDS-PFCO-1-EU</b>	ISOSAF Amber Barrel Enteral Syringe & Push Fit Cap, Low Dose, Single Use, BNS, 1ml	50pcs/ Pouch	500 units/box
<b>ABLDS-PFCO-3-EU</b>	ISOSAF Amber Barrel Enteral Syringe & Push Fit Cap, Low Dose, Single Use, BNS, 3ml	50pcs/ Pouch	500 units/box
<b>ABS-PFCO-5-EU</b>	ISOSAF Amber Barrel Enteral Syringe & Push Fit Cap, Single Use, BNS, 5ml	50pcs/ Pouch	500 units/box
<b>ABS-PFCO-10-EU</b>	ISOSAF Amber Barrel Enteral Syringe & Push Fit Cap, Single Use, BNS, 10ml	50pcs/ Pouch	500 units/box
<b>ABS-PFCO-20-EU</b>	ISOSAF Amber Barrel Enteral Syringe & Push Fit Cap, Single Use, BNS, 20ml	50pcs/ Pouch	200 units/box
<b>ABS-PFCO-35-EU</b>	ISOSAF Amber Barrel Enteral Syringe & Push Fit Cap, Single Use, BNS, 35ml	50pcs/ Pouch	100 units/box
<b>ABS-PFCO-60-EU</b>	ISOSAF Amber Barrel Enteral Syringe & Push Fit Cap, Single Use, BNS, 60ml	50pcs/ Pouch	100 units/box

# Feeding Tubes

## Nutricare Flo Feeding Tube

Nutricare Flo is an open ended polyurethane tube with clear numbered cm markings. It is fitted with the new ENFit (ISO 80369-3) connector.



### Product Features

- Rounded open end for enhanced flow rates
- 2 side port-holes
- 40cm - 125cm size range
- Radiopaque stripe
- Clearly numbered every CM
- Male ENFit connector
- Latex & DEHP free

ENFit®

Direct Codes	Product Descriptions	Packaging	Unit of Sale
NCF04.40	Nutricare Flo, No Stylet, Open End, Single Use, Nasogastric Feeding Tube, 4Fr x 40cm	1pc/pouch	25 units/box
NCF05.40	Nutricare Flo, No Stylet, Open End, Single Use, Nasogastric Feeding Tube, 5Fr x 40cm	1pc/pouch	25 units/box
NCF05.60	Nutricare Flo, No Stylet, Open End, Single Use, Nasogastric Feeding Tube, 5Fr x 60cm	1pc/pouch	25 units/box
NCF05.90	Nutricare Flo, No Stylet, Open End, Single Use, Nasogastric Feeding Tube, 5Fr x 90cm	1pc/pouch	25 units/box
NCF06.40	Nutricare Flo, No Stylet, Open End, Single Use, Nasogastric Feeding Tube, 6Fr x 40cm	1pc/pouch	25 units/box
NCF06.60	Nutricare Flo, No Stylet, Open End, Single Use, Nasogastric Feeding Tube, 6Fr x 60cm	1pc/pouch	25 units/box
NCF06.90	Nutricare Flo, No Stylet, Open End, Single Use, Nasogastric Feeding Tube, 6Fr x 90cm	1pc/pouch	25 units/box
NCF08.125	Nutricare Flo, No Stylet, Open End, Single Use, Nasogastric Feeding Tube, 8Fr x 125cm	1pc/pouch	25 units/box
NCF08.40	Nutricare Flo, No Stylet, Open End, Single Use, Nasogastric Feeding Tube, 8Fr x 40cm	1pc/pouch	25 units/box
NCF08.60	Nutricare Flo, No Stylet, Open End, Single Use, Nasogastric Feeding Tube, 8Fr x 60cm	1pc/pouch	25 units/box
NCF08.90	Nutricare Flo, No Stylet, Open End, Single Use, Nasogastric Feeding Tube, 8Fr x 90cm	1pc/pouch	25 units/box

\*Blue stripe on tube is for illustrative purposes only



# Enteral pH

## Enteral pH 2-9

Enteral pH 2-9 pH indicator strips.

### Product Features

- CE marked for the determination of human gastric aspirate
- 4 year shelf life
- Increments of 0.5
- Clear graduations between pH steps especially between 4.5-5.0 & 5.0-5.5
- pH range 2-9
- Can be stored from 40°C to 30°C

Available in packs of 25 for single patient use



### Patient Safety Alert

- In line with Patient Safety Alert NPSA/2011/PSA002 - Our pH indicator paper is CE marked and manufactured to test human gastric aspirate
- Supporting "Paediatric Nasogastric Tube Placement and Verification: Best Practice Recommendations From the NOVEL Project"

Irving, S. Y., Rempel, G. , Lyman, B. , Sevilla, W. M., Northington, L. , Guenter, P. and , (2018), Pediatric Nasogastric Tube Placement and Verification: Best Practice Recommendations From the NOVEL Project. Nutrition in Clinical Practice, 33: 921-927. doi:10.1002/ncp.10189

Direct Codes	Product Descriptions	Unit of Sale
ENT-PHSPPX	Enteral pH 2-9, Single Use, Tube of 25	25 strips/tube
ENT-PHX	Enteral pH 2-9, Single Use, Tube of 100	100 strips/tube



## pHX

pHX is a simple and effortless method for ascertaining the pH value of solutions with a color scale you are familiar with.

### Product Features

- CE marked for the determination of human gastric aspirate
- Increments of 0.5
- Accurate to pH +/- 0-2
- Supports Patient Safety Alert NPSA/2011/PSA002: Reducing the harm caused by misplaced nasogastric feeding tubes in adults, children and infants
- 0-6 & 2-9
- High quality pH paper
- 100 strips per container



Direct Codes	Product Descriptions	Unit of Sale
pHX0-6	Enteral pH 0-6, Single Use, Tube of 100	100 strips/tube
pHX2-9	Enteral pH 2-9, Single Use, Tube of 100	100 strips/tube

# ENPLUGs

## ENPLUG™ Stoma Plugs

ENPLUGs are designed for emergency use to keep stomas from closing when a gastrostomy device has fallen out.

The ENPLUG packs are available in 2 different lengths and various French sizes in case the stoma has already started to close and the current French gauge tube no longer fits the stoma.



### Product Features

- Stops stoma sites from closing
- Saves the stoma site until a new gastrostomy tube can be placed
- Includes a range of French gauge sizes in each pack which provides options even when stoma has started to close
- Effective, safe and inexpensive
- Saves staff time in emergency situations
- Promptly dealing with gastrostomy emergencies prevents use of other incorrect or inappropriate feeding devices
- Made from high grade medical silicone for flexibility
- Tapered tip for ease of insertion

### Pack Sizes

- **ENPLUG-S** Silicone Stoma Plug 4cm 10Fr, 12Fr, 14Fr, 16Fr
- **ENPLUG-L** Silicone Stoma Plug 7cm 10Fr, 12Fr, 14Fr & 16Fr
- **ENPLUG-XL** Silicone Stoma Plug 7cm 18 & 20Fr



Direct Codes	Specification	UOS
EN-PLUG-S	ENPLUG Silicone Stoma Plug Set, 4cm, 10Fr, 12Fr, 14Fr & 16Fr	100 units
EN-PLUG-L	ENPLUG Silicone Stoma Plug Set, 7cm, 10Fr, 12Fr, 14Fr & 16Fr	100 units
EN-PLUG-XL	ENPLUG Silicone Stoma Plug Set, 7cm, 18Fr & 20Fr	100 units



# Oral Syringes

## CareTip Oral Syringes

CareTip Oral Syringes available in single use for the delivery of feed/medications orally. Available in a variety of sizes with smooth barrel and tip.



1ml




2.5ml



5ml



10ml



20ml

### Product Features

- Single use
- Clear barrel markings
- 1ml with extra dead space reducing feature
- Single O ring for smooth plunging action
- Latex & DEHP free

Direct Codes	Specification	UOS
OS-1	CareTip™ Oral Syringe, Single Use, 1ml	100 units
OS-2	CareTip™ Oral Syringe, Single Use, 2.5ml	100 units
OS-5	CareTip™ Oral Syringe, Single Use, 5ml	100 units
OS-10	CareTip™ Oral Syringe, Single Use, 10ml	100 units
OS-20	CareTip™ Oral Syringe, Single Use, 20ml	100 units

# BNS Amber Oral Syringes

## Oral Syringes

Amber oral syringes with white markings for light sensitive medications. Amber barrel syringes do not transmit more than 10% of incident light with wavelengths of 290-450nm. As per USP671.



1ml



3ml



5ml



10ml



20ml

### Product Features

- Available as single use
- Clear barrel markings
- Latex & DEHP free

### With Push Fit Caps



Direct Codes	Specification	UOS
ABOSPFCO-1-500NS	Amber Barrel Oral Syringe & Push Fit Cap, Single Use, BNS, 1ml	500 units
ABOSPFCO-3-500NS	Amber Barrel Oral Syringe & Push Fit Cap, Single Use, BNS, 3ml	500 units
ABOSPFCO-5-500NS	Amber Barrel Oral Syringe & Push Fit Cap, Single Use, BNS, 5ml	500 units
ABOSPFCO-10-500NS	Amber Barrel Oral Syringe & Push Fit Cap, Single Use, BNS, 10ml	500 units
ABOSPFCO-20-200NS	Amber Barrel Oral Syringe & Push Fit Cap, Single Use, BNS, 20ml	200 units



# Bottle Adaptors

## ENFit Bottle Adaptor

Available in five sizes to fit the majority of medicine bottles.



### Product Features

- Available in 5 sizes to fit the majority of bottles
- Designed to fit tightly in the neck of the medicine bottle, allowing the bottle cap to be replaced
- Last the lifetime of the medicine bottle
- Purple (Enteral) color
- Latex & DEHP free
- Minimizes medicine waste
- Improve dose accuracy
- Available in single blisters or bulk packaging

Direct Codes	Specification	UOS
BAISO0-S	ISOSAF™ Medicine Bottle Adaptor & Tether, Single Use, Size 0: 8-10mm	100 units
BAISO1-S	ISOSAF™ Medicine Bottle Adaptor, Single Use, Size 1: 12-15mm	100 units
BAISO2-S	ISOSAF™ Medicine Bottle Adaptor, Single Use, Size 2: 14-16.5mm	100 units
BAISO3-S	ISOSAF™ Medicine Bottle Adaptor, Single Use, Size 3: 17-20.5mm	100 units
BAISO3.5-S	ISOSAF™ Medicine Bottle Adaptor, Single Use, Size 3.5: 19-23mm	100 units
BAISO4-S	ISOSAF™ Medicine Bottle Adaptor, Single Use, Size 4: 24-26mm	100 units
BAISO1	ISOSAF™ Medicine Bottle Adaptor, Single Use, Size 1: 12-15mm, Bulk	100 units
BAISO2	ISOSAF™ Medicine Bottle Adaptor, Single Use, Size 2: 14-16.5mm, Bulk	100 units
BAISO3	ISOSAF™ Medicine Bottle Adaptor, Single Use, Size 3: 17-20.5mm, Bulk	100 units
BAISO3.5	ISOSAF™ Medicine Bottle Adaptor, Single Use, Size 3.5: 19-23mm, Bulk	100 units
BAISO4	ISOSAF™ Medicine Bottle Adaptor, Single Use, Size 4: 24-26mm, Bulk	100 units

## Oral Bottle Adaptors

Bottle adaptor with oral connector.



### Product Features

- Available in 5 sizes to fit the majority of bottles
- Designed to fit tightly in the neck of the medicine bottle, allowing the bottle cap to be replaced
- Lasts the lifetime of the medicine bottle
- Latex & DEHP free
- Minimizes medicine waste
- Improves dose accuracy

Direct Codes	Specification	UOS
OBA-A	CareTip™ Bottle Adaptor for Oral Syringes, Single Use, Size A: 12-16mm	100 units
OBA-B	CareTip™ Bottle Adaptor for Oral Syringes, Single Use, Size B: 17-20mm	100 units
OBA-C	CareTip™ Bottle Adaptor for Oral Syringes, Single Use, Size C: 21-23mm	100 units
OBA-D	CareTip™ Bottle Adaptor for Oral Syringes, Single Use, Size D: 24-28mm	100 units



# Medicine Straws

## Medicine Straws

With Male ENFit connector (ISO 80369-3) available in PU or PVC material. Used to draw up feed/medication into an ENFit syringe. This allows for more accuracy, as per ISO80369-3 recommendations.

- Product Features
- Available in various lengths
  - Purple ENFit connector
  - Minimizes medicine waste
  - Improves dose accuracy
  - Longer straws avoid inverting larger medicine bottles
  - Latex & DEHP free
  - Ergonomic ridges on connector to allow easy connection/disconnection



Direct Codes	Specification	UOS
MSISO-50	ISOSAF™ Medicine Straw PVC, 6Fr x 50mm	100 units
MSISO-50P	ISOSAF™ Medicine Straw PU, 8Fr x 50mm	30 units
MSISO-100P	ISOSAF™ Medicine Straw PU, 8Fr x 100mm	30 units
MSISO-150P	ISOSAF™ Medicine Straw PU, 8Fr x 150mm	30 units

ENFit®



# Filter Needles

## Filter Needles

Used to draw up and filter medication into an ENFit Syringe via a glass ampoule or vial.

- Product Features
- Assembled with a protective sheath
  - 5µm filter for filling from glass ampoules
  - 18 gauge stainless steel needle with an angled tip allows for easy penetration through a vial's rubber membrane
  - 38mm & 50mm long for easy use with most vials
  - Commonly used with oral/enteral medications such as caffeine
  - Features an ENFit connection, allowing ENFit syringe users to draw up oral/enteral medications from manufacturers' vials, ampoules, and other containers
  - Ergonomic ridges on connector to allow easy connection/disconnection
  - Latex & DEHP free



Direct Codes	Specification	UOS
FNISO-38	ISOSAF™ Filter Needle, 5µm, 18G x 38mm	30 units
FNISO-50	ISOSAF™ Filter Needle, 5µm, 18G x 50mm	30 units

ENFit®





# Enteral Connector Brush

## Enteral Connector Brush

Enteral Connector Brush for keeping ENFit connectors free of debris.



### Product Features

- ENFit connectors have larger and deeper grooves, meaning increased area for debris to collect and potentially cause connection issues
- Using a brush after every use helps to prevent the accumulation of debris, and the feeding tube cap from adhering to the connector
- Keeping the connector and cap free of debris can help to extend the lifespan of the tube

In line with GEDSA ENFit cleaning Procedures 2018 recommendations Ref <http://stayconnected.org/wp-content/uploads/2018/05/ENFit-Cleaning-Procedures-All-Tubes.pdf>



Direct Codes	Specification	UOS
ECB-NS	Enteral Connector Brush	50 units




# Clear Grip

## Clear Grip

Clear Grip unique skin fixation system designed specifically to hold tubing in place. E.g. NG Feeding extensions.

### Product Features

- Secures tubing between 2 layers of material, away from the patient's skin
- Clear material allows observation of skin
- Ergonomic design ensures a comfortable fit, close to nostril to minimise tube pull outs
- For use in pediatric patients and older use, not for NICU patients



Direct Codes	Specification	UOS
CG-30	Clear Grip Fixation Plaster	30 units
CG-NS	Clear Grip Fixation Plaster, Bulk	2000 units



# Enteral Caps

## Enteral Caps

Enteral Cap Male & Female  
ENFit connector (ISO 80369-3).



ENFit Syringe  
Cap



ENFit Female  
Cap



ENFit Female  
Cap

### Product Features

- Secure syringe caps to allow medication to be prepared in advance
- Female ENFit caps close off extension tubing or similar female ENFit connectors

Direct Codes	Specification	UOS
ECISO	ISOSAF™ Male Syringe Cap, Single Use	100 units
FCISO-10	ISOSAF™ Female Cap, Single Use	150 units
FCISO-NS	ISOSAF™ Female Cap, Single Use	200 units

ENFit®



# Push Fit Caps

## Push Fit Caps

Male cap for both  
enteral and oral syringes.



### Product Features

- Auto stand cap
- Syringe cap allows medication to be prepared in advance
- Push fit mechanism for ease of use
- Allows syringes to be capped off until required
- Self righting cap, allows for one handed capping

Direct Codes	Specification	UOS
PFC-100	Push Fit Cap - Enteral & Oral, Single Use	100 units
PFC-500NS	Push Fit Cap - Enteral & Oral, Single Use, BNS	500 units
PFCO-200NS	Orange Push Fit Syringe Cap, Single Use, BNS	200 units
PFCO-500NS	Orange Push Fit Syringe Cap, Single Use, BNS	500 units

## Stepped Bottle Adaptor

Stepped Bottle Adaptor,  
16-28mm, male ENFit  
connector (ISO 80369-3).



### Product Features

- Universally sized bottle adaptor
- Features both ENFit and oral connector

Direct Codes	Specification	UOS
SBAISO	ISOSAF™ Stepped Medicine Bottle Adaptor, Single Use	50 units




# Drainage Bags

## Enteral Drainage Bags


Enteral Drainage Bags to fit any ENFit compliant tubes, are used for gastric drainage.

The bags are available in either soft EVA material (50ml - 1000ml) or PVC Material (250ml - 2000ml)



### Product Features

- Clamp
- Latex & DEHP free
- Tethered cap
- Male ENFit & funnel connector available
- ENFit compliant



ENFit®

	Direct Codes	Specification	UOS
EVA Drainage Bags	DBISO-50	EVA Drainage Bag, Single Use, 50ml	20 units
	DBISO-100	EVA Drainage Bag, Single Use, 100ml	20 units
	DBISO-250	EVA Drainage Bag, Single Use, 250ml	20 units
	DBISON-500	EVA Drainage Bag, Single Use, 500ml	20 units
	DBISO-500	EVA Drainage Bag with Funnel Connector, Single Use, 500ml	20 units
	DBISO-1000	EVA Drainage Bag with Funnel Connector, Single Use, 1000ml	20 units

	Direct Codes	Specification	UOS
PVC Drainage Bags	DBP-250	PVC Drainage Bag, Single Use, 250ml	50 units
	DBP-500	PVC Drainage Bag, Single Use, 500ml	50 units
	DBP-1000	PVC Drainage Bag, Single Use, 1000ml	50 units
	DBP-1500	PVC Drainage Bag, Single Use, 1500ml	50 units
	DBP-2000	PVC Drainage Bag, Single Use, 2000ml	50 units

## Notes

[illegible]This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.





GBUK Group award-winning specialist companies are experts in the provision of enteral feeding, critical care and patient moving and handling products and technologies. GBUK Group companies are suppliers to every UK NHS Hospital and to a rapidly expanding export customer base.

**GBUK Group Ltd.**

Woodland House, Blackwood Hall Business Park,  
North Duffield, Selby, North Yorkshire, YO8 5DD, UK

T +44 (0) 1757 288 587

E [info@gbukgroup.com](mailto:info@gbukgroup.com)

W [www.gbukgroup.com](http://www.gbukgroup.com)

ENT00015\_1.0 2012 05