




By Pass Adapter

To maintain a closed circuit



Integrated Swing Valve Technology™



Made in Canada

Flusso™ and Swing Valve Technology™ are trademarks used under licence by
McArthur Medical Sales Inc.

2018

SILVER WINNER

**MEDICAL DESIGN
EXCELLENCE AWARDS**

US and International Patents Pending

Flusso™ Highlights:

- Proprietary valve with Swing Valve Technology™
- Clear housing with colored valve permits visualization of function of Swing Valve Technology™
- Integrated swivel at patient connection to reduce torque on endotracheal tube
- Tethered dust cap keeps ports clean and protects swing valve
- Minimal dead space (9cc)
- Standard ISO 15/22mm connections

Flusso facilitates the disconnection of the patient from the mechanical ventilator during a circuit change, HME change, or disconnection of patient for transport.

Flusso By Pass allows for the maintenance of a patient's positive end expiratory pressure and reduce potential exhaled aerosols during ventilator circuit disconnect.

Product Code: 420-1978; 25/case



www.keepthepeep.com

A Product from:



Solutions in Critical Care

Manufactured by:



Distributed by:



For More Information and a Sample

www.thinkmdi.net

866-624-3952