




To Enable a Safer Circuit Disconnect



Made in Canada



MMSI
McARTHUR
MEDICAL SALES INC.



Flusso™ TFI allows for temporary interruption of gas flow to maintain a patient's positive end expiratory pressure (PEEP) and minimize potential exhaled aerosols during ventilator circuit disconnect.



www.keepthepeep.com

FLUSSO™
TFI

Flusso™ Highlights:

- Lightweight
- Clear housing permits visualization of function
- Low deadspace (1.86 cc)
- Compact
- Minimal effort to activate
- Safety tabs to prevent over compression

Product Code: 420-1979; 25/case

Manufactured by:



A Product from:



Solutions in Critical Care

Distributed by:

Medical Dynamics
TECHNOLOGY OF THE FUTURE, TODAY

For More Information and a Sample
www.thinkmdi.net • 866-624-3952

Premier Contract Number: PP-DS-743
Vizient Contract Number: DM0170