

## Right-Fit Pediatric Patients with Excellent Seal and Comfort

The Artemis Pedi-Fit Non-Invasive Pediatric Nasal and Full-Face Masks were specifically designed to right-fit respiratory patients in the Pediatric Intensive Care Unit (PICU). The soft, molded silicone material and flexible forehead pad ensure excellent seal and patient comfort during CPAP and VPAP therapy.

### Pedi-Fit Features

- Molded silicone mask cushion provides excellent seal and comfort
- Unique clip pivots the mask to conform to the patient's nose, preventing air from leaking into the patient's eyes
- Flexible, one-piece forehead pad evenly conforms to the forehead for stability
- Quiet air vent minimizes noises and enables efficient CO<sub>2</sub> washout
- Each design is available in four sizes enabling healthcare providers to right-fit each patient

Item ID	Description	U/M
ART00PF	Pedi-Fit Nasal Mask, Small (21 mm x 18.6 mm)	Each
ART0PF	Pedi-Fit Nasal Mask, Medium (22 mm x 24.10 mm)	Each
ART01PF	Pedi-Fit Nasal Mask, Large (26.5 mm x 28.2 mm)	Each
ART02PF	Pedi-Fit Nasal Mask, Adult X-Small (30.5 mm x 32.3 mm)	Each
ART200PF	Pedi-Fit Full-Face Mask, Small (28mm x 37 mm)	Each
ART20PF	Pedi-Fit Full-Face Mask, Medium (35 mm x 43 mm)	Each
ART21PF	Pedi-Fit Full-Face Mask, Large (41 mm x 50 mm)	Each
ART22PF	Pedi-Fit Full-Face Mask, X-Large (46 mm x 60 mm)	Each
TST600	72" Ventilator Circuit with Smooth Bore Tubing	20/Bx



**Nasal Mask**



**Full-Face Mask**



**Ventilator Circuit**

Distributed By

**Medical Dynamics**  
TECHNOLOGY OF THE FUTURE, TODAY


5720 Lonetree Blvd., Rocklin, CA 95765  
www.thinkmdi.net • (866) 624-3952




**Solutions in Critical Care**

Premier Contract Number: PP-DS-743  
Vizient Contract Number: DM0170

### Sizing Chart for Pedi-Fit Nasal Mask



### Sizing Chart for Pedi-Fit Full-Face Mask



Distributed By

**Medical Dynamics**  
TECHNOLOGY OF THE FUTURE, TODAY

5720 Lonetree Blvd., Rocklin, CA 95765  
www.thinkmdi.net • (866) 624-3952



**Solutions in Critical Care**

Premier Contract Number: PP-DS-743  
Vizient Contract Number: DM0170