

# Sanitize Staff & Visitor Devices At The Speed of Light.



An average of **1 in 4 hospital mobile devices are contaminated with pathogenic bacteria.**

**We can help fix this.**



## A Safe & Easy Solution to Reduce the Spread of Bacteria on Mobile Devices.



Fast and simple.  
No training required.



Proven to inactivate  
99.9998% of MRSA  
without pre-cleaning.



Touch-free device  
removal that enables  
proper hand hygiene.



No chemicals and no  
damage to sensitive  
touchscreens.

**Break the Chain of Infection**

**Designed for Healthcare.**  
Simple, fast, and effective  
device sanitization.



**Designed for Hospitals.  
Trusted Everywhere.**

**2000+**  
Deployments

**Zero**  
Damage to  
Touchscreens

**300+**  
Hospitals Across  
USA & CAN

**720**  
Devices per  
Hour Capacity.



### Pathogenic Kill Rates @ 20 Secs

- *S. aureus* (MRSA) - 99.99981% reduction
- *C. difficile* spores - 99.97% reduction
- *S. enterica* (Salmonella) - 99.999% reduction
- *E. faecalis* (VRE) - 99.9992% reduction
- *E. coli* - 99.9991% reduction
- SARS-CoV-2 - 99.995%

### Specifications

- **Size:** 26.7" x 16.5" x 12.3"
- **Tray:** 9.4" x 12.87" x 3.45"
- **Desktop Footprint:** 15.7" x 16.45"

**Trusted by Over  
300 Hospitals**  
to Sanitize Phones,  
Tablets & More.

"The CleanSlate is getting a workout in our unit... it has truly become integrated into our patient safety protocols. The families and staff LOVE it. Anything that helps keep our tiny miracles fighting another day, is such a blessing."  
- Stephanie P., NICU Manager, UMPC Pinnacle

**cleanslate** UV

Provided By:  
**Medical Dynamics**  
TECHNOLOGY OF THE FUTURE, TODAY™  
866-624-3952  
www.thinkMDI.net

**USA**  
1170 Main St  
Buffalo, NY USA  
14209

**sales@cleanslateuv.com**