



# Sanitize Staff & Visitor Devices. Protect NICU Patients.

Staff and visitors are required to wash their hands as they enter Neonatal Intensive Care (NICU) - but what about the devices in their pockets, around their necks or on their belts?

These devices are among the dirtiest surfaces in the hospital, and are often touched immediately after we wash or sanitize our hands.

CleanSlateUV's 20 second, chemical-free solution allows NICUs to address this problem. The intuitive solution can be placed next to hand washing stations in your NICU, so staff and visitors can sanitize their hands and devices with ease.

# Designed for Healthcare.

Simple, fast, and effective device sanitization.



## The Trusted Solution.

Used in 300+ Hospitals with 50+ NICU Deployments.

Five "U.S. News: Best Children's Hospitals Honor Roll" Clients

4,000,000+ NICU devices sanitized.

### A simple 20-second process

proven to inactivate 99.9998% of MRSA without pre-cleaning.

### Prevent re-contamination

of clean hands with touch-free device removal.

### Works on a range of devices,

including smartphones, tablets, voceras, badges, and other non-critical devices.

### No device-damaging chemicals

interacting with touchscreen devices or expensive equipment.

## Pathogenic Kill Rates @ 20 Secs

- *S. aureus* (MRSA) - 99.99981% reduction
- *C. difficile* spores - 99.97% reduction
- *S. enterica* (Salmonella) - 99.999% reduction
- *E. faecalis* (VRE) - 99.9992% reduction
- *E. coli* - 99.9991% reduction
- SARS-CoV-2 - 99.995%

## Specifications

- **Size:** 26.7" x 16.5" x 12.3"
- **Tray:** 9.4" x 12.87" x 3.45"
- **Desktop Footprint:** 15.7" x 16.45"

## Take the First Step to Help Your NICU

"The CleanSlate is getting a workout in our unit... it has truly become integrated into our patient safety protocols. The families and staff LOVE it. Anything that helps keep our tiny miracles fighting another day, is such a blessing."

- Stephanie P., NICU Manager, UMPC Pinnacle

"Every time we bring a patient into the NICU, we need to tell their parents not to use their phones. We don't have any way to make sure the phones won't become a vector for transmission. And we usually have to tell the parents this several times."

- Jean C., Catholic Health NICU

**cleanslate** UV

Provided By:  
**Medical Dynamics**  
TECHNOLOGY OF THE FUTURE, TODAY™  
866-624-3952  
www.thinkMDI.net

**USA**  
1170 Main St Buffalo, NY, USA  
14209