

# DimAir Line



## TOTAL FACE DIMAX

Total Face mask for non-invasive mechanical ventilation in PSV, BIPAP and CPAP techniques. Featuring 5 attachment points, it optimizes patient/ventilator synchronization and pressure distribution to the face, reducing leakages and improving comfort. The high volume dedicated to the nose and the possibility of replacing the connectors to transform the same mask into Vented or Non Vented, make DIMAX the ideal interface for non-invasive ventilation.

**7 plus**

- + Connection
- + Comfort
- + Sealing
- + Practicality
- + Flexibility
- + Safety
- + Versions

## TECHNICAL FEATURES:

- 5 fixation points.
- More space for the nose.
- Dedicated access for NasoGastric tubes
- Interchangeable Connectors:  
22M Vented  
22F Non Vented

## ADVANTAGES:

- Better patient/ventilator synchronism.
- Improved pneumatic seal.
- More comfort for prominent noses.
- Lower headgear traction.
- Reduced risk of contact injury.



The 5th fastening point allows greater grip on the patient's forehead



Kit with Easy Vent



"Non Vented" connector with anti asphyxia valve



Dedicated passage for double-lumen tube



"Non Vented" Connector for bitube ventilator



"Vented" connector with anti asphyxia valve for single tube ventilator



## REFERENCES:

Item Code	Description	U/M
700/7309US	DIMAX Total Face Mask, Small, headgear, ngt port, three elbows	5/box
700/7310US	DIMAX Total Face Mask, Medium, headgear, ngt port, three elbows	5/box
700/7311US	DIMAX Total Face Mask, Large, headgear, ngt port, three elbows	5/box
700/7312US	DIMAX Total Face Mask, XL, headgear, ngt port, three elbows	5/box
700/7313US	DIMAX Total Face Mask, XXL, headgear, ngt port, three elbows	5/box



**HMC Premedical S.p.A**

Via Bosco 1/3 – 41037 Mirandola – MO – I  
Tel: +39 0535 22704 Fax: +39 0535 609546



**Solutions in Critical Care**

www.solutionsincriticalcare.com  
866-367-7221